TECHNIQUE CHIRURGICALE

HERMES Revision

LE GENOU DE REVISION PEU CONTRAINT

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HERMES Revision

THE LOW STRESS REVISION KNEE

- Easy to use instrumentation to face all situations
- Stability of cuts
- Preservation of bone
- Re-alignment of the mechanical axis in extension
- Ligament balance in extension and in flexion

SURGICAL TECHNIQUE: STANDARD TIBIAL CUT AND TIBIAL TRIAL IMPLANT

Tray #50281

1 ASYMETRIC TIBIAL CUT

Aim :
Asymetric tibial implant (spacer of 10 mm)

2 TRIAL & TIBIAL BASE PLATE

Tibial implant with a 10 mm spacer

Asymetric tibial implant (spacer of 10 mm) internal or external side

3 TRIAL & TIBIAL BASE PLATE HIFIT

HiFit Standard tibial implant

		HiFit PS FEMORAL COMPONENT								
		1	2	3	4	5	6	7	8	9
	Femoral Component HERMES PS REVISION	1	2	3	4	5	6	7	8	9
Tibial component	1	✓	✓	✓	✓	✓				
	2	✓	✓	✓	✓	✓				
	Peg 3	✓	✓	✓	✓	✓				
	Peg 2						✓	✓	✓	✓
	Peg 3						✓	✓	✓	✓
	4						✓	✓	✓	✓

		FEMORAL COMPONENT HERMES PS REVISION			
		1	2	3	4
HiFit PS TIBIAL INSERT	A	✓ Small peg	✓ Small peg		
	B	✓ Small peg	✓ Small peg		
	C	✓ Small peg	✓ Small peg	✓ Large peg	✓ Large peg
	D	✓ Small peg	✓ Small peg	✓ Large peg	✓ Large peg
	E			✓ Large peg	✓ Large peg
	F			✓ Large peg	✓ Large peg

SURGICAL TECHNIQUE: AP FEMORAL CUTS AND TRIAL FEMORAL COMPONENT

Tray #50282

4 ANTERO-POSTERIOR FEMORAL CUTS

Aim :
• Femoral stem diameter and angle (10 to 20 mm ; 5° or 7°)
• Femoral component size (1 to 4)
• Posterior femoral spacer (5 or 10 mm)

5 DISTAL FEMORAL CUT

Aim :
Distal femoral spacer (5 or 10 mm)

6 MEASUREMENT OF FEMUR TIBIA GAP

Aim :
Tibial inlay thickness (7, 9, 11, 13, 15, 17)

Compensate the asymetric tibial cut with a 10 mm spacer

7 PREPARING THE FEMUR

8 TRIALS

9 ASSEMBLING THE EXTENSION OR ADAPTOR TIBIAL EXTENSION

Unscrewing the screw in the base to release the PE plug

Screwing the tibial extension or eccentric adaptor with the tibial extension until tightening

Insert three layers of compresses, then hammer three times while holding the base plate on the table

Finish screwing the screw into the extension or eccentric

PE plug

Tibial extension or eccentric adaptor with the tibial extension

Dimensional correspondence of femoral components HiFit and HERMES PS revision

Femoral HiFit	1	2	3	4	5	6	7	8	9
Femoral HERMES PS Revision	1			2				4	

CERAVER - LES LABORATOIRES OSTÉAL MEDICAL
69, RUE DE LA BELLE ÉTOILE - 95957 ROISSY CDG CEDEX - FRANCE
TEL. : (+33) 1 48 63 88 63 - FAX : (+33) 1 49 89 05 05 - WWW.CERAVER.COM - CONTACT@CERAVER.COM

TOG3703.20