● NOTICE D'INSTRUCTIONS

Français

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● INSTRUCTIONS FOR USE

English

Instructions for use for the following implants:

- Stainless steel heads (References: 2220, 2221, 2222, 2223, 2263, 2264, 2223, 2224, 2225, 2261, 2262, 2217, 2218, 2219)

Graphic symbols used for labelling these devices:

Symbol	Symbol title	Symbol	Symbol title
	Manufacturer		Do not reuse
	Expiry date		Product terminally sterilized by irradiation.
	Batch code		Indicates a medical product having undergone sterilisation.
	Catalogue reference of the device		Caution, all the information required for the use of the device is not found on the label; it is necessary to consult the instructions for use.

Before using a CERAVER® hip prosthesis, the attention of the surgeon is drawn to the following recommendations:

1 - GENERAL POINTS

- 1.1 The manufacturer and supplier of a CERAVER® femoral head guarantee the quality of the manufacture and materials. A CERAVER® femoral head that has already been used must not be re-implanted. Instructions for the assembly of the different components must be adhered to as defined in the commercial and technical operating documentation. The use of CERAVER® prosthetic hip components with components of other origins is strictly prohibited.
- 1.2 Several criteria are of crucial importance for the success of an operation using a CERAVER® femoral head:
 - Appropriate patient selection for the relevant procedure
 - Perfect operating technique and correct use of the surgical instruments specially designed by CERAVER®.
 - The strictest asepsis.
- 1.3 The surgeon is responsible for any complications that may arise from an erroneous indication and/or a defective operating technique or incorrect asepsis; these cannot be attributed to the manufacturer or the supplier of the CERAVER® products.
- 1.4 The surgeon must be aware of all problems posed by implanting a joint prosthesis. The operating technique for implantation of a CERAVER® hip prosthesis may be taught in a surgical centre where surgeons have perfect mastery of this technique.
- 1.5 Patients receiving hip prostheses should be informed that their weight and activity level may influence the durability of the prosthesis.

2 - IMPLANT DESCRIPTION

CERAVER® femoral heads are manufactured using the latest technologies. The mechanical properties and functional characteristics can only be guaranteed if the surgical technique and present instructions are strictly observed.

- The femoral heads are made of stainless steel FeCrNiMnMo with high nitrogen content (ISO 5832-9 standard, ASTM F 1586 standard)

3 - INDICATIONS

- The implantation of a hip prosthesis may be considered when all other surgical possibilities have been carefully weighed up and are considered to be less appropriate.
- Even the most successfully implanted artificial joint is inferior to a sound natural joint. On the other hand, an artificial joint can be a highly beneficial substitute for a deformed or diseased joint, because it eliminates pain and is conducive to the restoration of good mobility and weight-bearing capacity of the operated limb.
- Every artificial joint is always subject to wear which cannot be avoided. An initially stable prosthesis may become loose in the course of time. Wear and loosening of a prosthesis may render a re-operation necessary.
- Infection of a prosthesis will usually necessitate its removal and any such infection may leave debilitating sequelae.

Approved applications:

Hip arthroplasty involves the total or partial replacement of a hip joint affected by a variety of conditions, which cannot be treated by conservative treatment or another surgical treatment. Its purpose is to give the patient back all or part of his/her mobility and stability while minimising or even eliminating the pain.

A hip joint replacement may be indicated in the following cases:

- primary hip osteoarthritis
- osteonecrosis of the femoral head
- trauma sequelae
- epiphysiolysis
- hip dysplasia
- hip dislocation
- sequelae of hip infections
- fracture of the proximal femur
- non-infectious arthritis
- revision surgery (arthrodeseis) or total hip replacement.

In general, patients must be suffering from severe joint involvement which is unlikely to respond to conservative treatment or another surgical strategy.

Given the considerations mentioned above, hip replacement is advised when the following conditions are met:

- The anatomical shapes and quality of the bone structures in patients must allow the anchorage of a prosthesis.
- The surgeon must inform patients of the risks of implanting a prosthesis and their consent to the proposed operation must be given.

Remark: In the event of revision surgery after fracture of the alumina head, given the presence of particles in the peri-prosthetic tissues, as far as possible it is preferable to re-implant an alumina head rather than a metal head, in order to prevent any risk of abrasion of the metal head.

The following circumstances require special precautions:

- Overweight patient.
- Excessive weight-bearing expected due to his/her profession or participation in sport.
- Epilepsy or other factors that may generate repetitive accidents with an increased risk of fractures.
- Severe osteoporosis and/or osteomalacia.
- Potential infectious diseases with potential for evolution, with possibility of joint symptoms.
- Significant deformation of the joint to be replaced, indicating possible difficulties in anchoring the prosthesis.
- Reduction in bone support capacity following a tumour.
- Intellectual inability of patients to understand and follow the surgeon's instructions.
- Alcoholism or drug addiction.
- Pre-existing infectious or general drug treatments leading to progressive deterioration of the bone support for the prosthesis.
- Disorders in other joints (knee, ankle, etc.).
- Globally weakened immune status of the patient.
- Unbalanced diet.

4 - CONTRAINDICATIONS

The following circumstances generally constitute an absolute contraindication to the implantation of a joint prosthesis:

- Allergy to any component of the prostheses.
- Active infection of the joint or relevant peri-articular area
- Loss of bone substances above or below the joint, making it impossible to anchor the cementless prosthesis correctly.
- Charcot disease, severe neuroarthropathies, muscle loss, neuromuscular lesion or vascular insufficiency of the affected limb, making the operation unjustifiable.
- Patients under 50 years old suffering from joint lesions that could probably be treated successfully by a conservative treatment.
- Pregnancy
- Do not combine a ø28 or ø32 stainless steel femoral head, neck +12 (references 2262 and 2264) with a CERAFIT R-MIS lateralized stem size 16 (references 12416 and 12216). This configuration does not meet the requirements of the current NF ISO 7206-6 standard (the phenomenon observed is a fatigue failure of the neck in less than 10 million cycles). This configuration is therefore prohibited under penalty of exposing the patient to a proven risk of fatigue failure of the neck of the femoral stem.

5 - IMPORTANT WARNINGS

5.1 Preoperative warnings

- The operation must be scheduled with care. Depending on the results of the X-ray examination (templates) and preoperative planning, the surgeon will decide on the suitable implant.
- CERAVER® instruments and templates specially designed to prepare for the implantation of the joint components should be used. The description of the surgical procedure can never be exhaustive or specify all the risks and complications that must be accounted for. Technical notes about our products and the appearance of functional manifestations. Early recognition of a complication enables effective measures to be taken. In particular, early re-operation has more chance of success than a late operation, particularly when advanced loosening has already led to significant reabsorption lesions in the bone.
- Patients under 50 years old suffering from joint lesions that could probably be treated successfully by a conservative treatment.
- Pregnancy
- Do not combine a ø28 or ø32 stainless steel femoral head, neck +12 (references 2262 and 2264) with a CERAFIT R-MIS lateralized stem size 16 (references 12416 and 12216). This configuration does not meet the requirements of the current NF ISO 7206-6 standard (the phenomenon observed is a fatigue failure of the neck in less than 10 million cycles). This configuration is therefore prohibited under penalty of exposing the patient to a proven risk of fatigue failure of the neck of the femoral stem.

5.2 Postoperative warnings

- The joint prostheses must be handled with care. Scratches or unevenness due to impacts will lead to wear and weight bearing capacity will be reduced.
- For the choice of prosthesis size: the templates, ancillary instruments and CERAVER® test implants will be the only ones used.
- Bare prosthetic joint surfaces must not be in contact with any metal object, instruments for example.

6 - POSTOPERATIVE PRECAUTIONS

6.1 All patients implanted with a CERAVER® joint prosthesis should see an orthopaedic surgeon for postoperative follow-up. This examination is important to detect any signs of prosthetic wear or loosening before the appearance of functional manifestations. Early recognition of a complication enables effective measures to be taken. In particular, early re-operation has more chance of success than a late operation, particularly when advanced loosening has already led to significant reabsorption lesions in the bone.

6.2 The patient must inform his/her surgeon about the slightest change occurring in the operated limb.

6.3 If an annual follow-up appointment is not possible, a control X-ray image should be sent to the surgeon for evaluation at least every year. The radiological study makes it possible to screen for possible early lesions.

6.4 MRI compatibility: The CERAVER® stainless steel heads are non-ferromagnetic, have a geometry that is unlikely to generate induced currents, and are fixed to the implant, which reduces the risk of mobilisation. Their safety, particularly in terms of implant heating and migration, was assessed from bibliographical data by comparison with the available data on devices of similar composition, shape, and use. This evaluation has concluded that there is no significant compatibility for MRI examinations up to 3 Tesla. However, as a precaution, it is recommended that MRI examinations are avoided in the 48 hours after implantation, and that the person responsible for the examination of the recent insertion of an implant is warned, if such an examination proves to be essential.

It should be noted that all devices presenting a high contrast with the surrounding biological environment may generate artefacts which must be taken into account for the perfect performance and interpretation of imaging examinations. For this purpose, a patient with an implant should be recommended to inform the health professionals concerned (radiologist and radiology technicians) about the presence of an implant as soon as possible before these examinations.

7 - COMPLICATIONS

The following complications have occasionally been observed, and require special attention:

- Loosening of the prosthesis following changes in the transmission of weight bearing or tissue reactions due to the prosthesis. Loosening is frequently the result of one or more of the risk factors previously mentioned.
- Early and late-onset infections.
- Dislocation, subluxation, insufficient range of motion, undesirable shortening or extension of the operated limb following a non-optimal implantation of the prosthesis.
- Too great a medialisation of the joint which could cause pain in other joints or the back.
- Bone fractures following unilateral overload or bone fragility.
- Transient or permanent nerve lesions following compression or haematoma.
- Cardiovascular diseases or disorders, including venous thromboses, pulmonary embolisms and cardiac arrhythmias. Damage to the blood vessels during the operation.
- Haematoma of the surgical wound and delay in healing.
- Tissue reactions following an allergy to the implanted materials, particularly metal, or following the accumulation of wear particles.
- Temporary or permanent peripheral neuropathy.
- Pulmonary disorders, particularly pneumonia and atelectasis.
- Urological complications (retention, urinary infection).
- Other associated complications resulting from anaesthesia, use of ancillary equipment, nosocomial diseases.
- Ectopic ossification.
- Pain.

8 - COMPATIBLE MEDICAL DEVICES

The stainless steel femoral heads are compatible with:

- The standard polyethylene liners of the same joint diameter.
- The PE cemented acetabular cups of the same joint diameter
- The PE dual mobility cups (only with femoral heads of Ø 28 mm)
- The bipolar heads of the same joint diameter (Ø 22,2 or 28 mm)
- All cementless stems in the CERAFIT range with the exception of the combination mentioned in paragraph 4 - Contraindications -.
- All cementless stems in the OSTEAL range
- Dedicated ancillary instrumentation

Comment: joint diameter means the nominal diameter of the bearing.

9 - STORING AND HANDLING PROSTHESES

- Joint prostheses must be stored in the original packaging.
- Joint prostheses must not be subjected to treatment or modification.
- Prostheses that are damaged or have undergone a treatment that is inappropriate or not authorized by the manufacturer must not be implanted.
- Do not use the device after the expiry date. Contact your CERAVER® representative and return the product to him/her.

10 - PACKAGING AND STERILISATION

The medical devices are packed individually in double packaging and are supplied sterile. CERAVER® femoral heads are sterilized by gamma radiation at a minimum dose of 25 kGy. It is important to check the integrity of the sterility protection before use.

11 - RESTERILISATION AND REUSE

Resterilisation and/or reuse of any of the components of a CERAVER® hip prosthesis is strictly prohibited at the risk of exposing the patient concerned to contamination that could cause an infection.

12 - DISPOSAL AFTER USE

If a device is removed for revision surgery, dispose of the components according to the centre's medical waste disposal procedure. This does not apply if the follow-up procedure is due to a Medical Device Vigilance case.

13 - NOTE

The stainless steel femoral heads were CE marked in 1998. Federal law in the USA and Canada requires that prosthetic components be purchased by a pharmacist. For further information concerning the use of these components, please contact the manufacturer, its overseas subsidiaries and its distributors.

CERAVER - Les Laboratoires Osteal Médical
69 rue de la Belle Etoile
95057 Roissy CDG Cedex - France
Tél.: 33 (0)1 48 63 88 63 - Fax: 33 (0)1 48 63 88 99
www.ceraver.com

Instructions last updated: February 2025

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