## ● NOTICE D'INSTRUCTIONS COTYLE

### Français

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## ● INSTRUCTIONS FOR USE OF CUP

### English

#### Instructions for use for the following implants:

- HAC cementless dual mobility cup with disposable gripper (ref.: 12019, 12020, 12021, 12022, 12023, 12024, 12025, 12026, 12027)

Graphic symbols used for labelling these products:

| Symbol                                                                            | Symbol title                      | Symbol                                                                             | Symbol title                                                                                                                                    |
|-----------------------------------------------------------------------------------|-----------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Manufacturer                      |  | Do not reuse                                                                                                                                    |
|  | Expiry date                       |  | Product terminally sterilized by irradiation.                                                                                                   |
|  | Batch code                        |  | Indicates a medical product having undergone sterilisation.                                                                                     |
|  | Catalogue reference of the device |  | Caution, all the information required for the use of the device is not found on the label; it is necessary to consult the instructions for use. |

Before using a CERAVER® hip prosthesis, the attention of the surgeon is drawn to the following recommendations:

#### 1 - GENERAL POINTS

1.1 The manufacturer and supplier of a CERAVER® hip prosthesis guarantee the quality of the manufacturing process and materials used. A prosthesis that has already been used must not be implanted again. Instructions for the assembly of the different components have to be respected, as specified in the commercial documentation. The use of CERAVER® prosthetic hip components with components of other origins is strictly prohibited.

1.2 Several criteria are key to the success of a procedure to implant a CERAVER® prosthesis:

- Appropriate patient selection for the relevant procedure.
- Perfect surgical technique and correct use of the surgical instruments specially designed by CERAVER®.
- Application of strict aseptic technique.

1.3 The surgeon is responsible for any complications that may arise from an erroneous indication and defective operating technique or incorrect asepsis; these cannot be attributed to the manufacturer or the supplier of the CERAVER® products.

1.4 The surgeon must be aware of all problems posed by implanting a joint prosthesis. The surgical technique for implantation of a CERAVER® hip prosthesis may be taught in a surgical centre where surgeons have perfect mastery of this technique.

1.5 Patients receiving a prosthetic hip joint should be informed that their weight and activity level may affect the longevity of the prosthesis.

#### 2 - IMPLANT DESCRIPTION

The components of CERAVER® cementless hip prostheses are manufactured using the latest technologies. The mechanical properties and functional properties can only be guaranteed if the surgical technique and present instructions are strictly observed. The HAC cementless dual mobility cup is designed to be implanted without cement.

- The HAC cementless dual mobility cup is made of stainless steel with high nitrogen content (ISO 5832-9 standard, ASTM F 1586 standard) coated with hydroxyapatite (HA) (ISO 13779-2, ISO 13779-3, ISO 13779-4 standards) on the side in contact with the acetabulum.

- The disposable gripper is made of ultra high molecular weight polyethylene (UHMWPE).

The range of angular motion permitted by CERAVER® dual mobility cups complies with the requirements of standard NF EN ISO 21535.

#### 3 - INDICATIONS

- A hip prosthesis should be considered only after all other surgical possibilities have been carefully considered and deemed to be less appropriate.
- Even the most successfully implanted artificial joint is inferior to a sound natural joint. On the other hand, an artificial joint can be a highly beneficial substitute for a deformed or diseased joint, because it eliminates pain and is conducive to the restoration of good mobility and weight-bearing capacity of the operated limb.
- Every artificial joint is always subject to wear which cannot be avoided. An initially stable prosthesis may become loose in the course of time. Wear and loosening of a prosthesis may render a re-operation necessary.
- Significant deformation of the joint to be replaced, indicating possible difficulties in anchoring the prosthesis.
- Infection of a prosthesis will usually necessitate its removal or replacement under certain conditions and any such infection may leave debilitating sequelae.

#### Approved applications:

Hip arthroplasty involves replacing a hip joint subject to a variety of conditions that cannot be treated by other therapies or methods. Its purpose is to give the patient back all or part of his/her mobility and stability while minimising or even eliminating the pain.

The clinical situations for use of the CERAVER® HAC cementless dual mobility cup system are as follows

- Degradation of the joint as a result of primary or secondary osteoarthritis (including avascular necrosis, or congenital hip dysplasia) or inflammatory coxopathy
- Significant fracture of the proximal end of the femur;
- Condition after previous surgical procedures such as joint reconstruction (osteotomy), arthrodesis or hemiarthroplasty.

For HAC cementless dual mobility cups, the specific indications concern:

- Risk of luxation associated particularly with a neuromuscular deficit, a dysplastic problem, or physical or mental impairment of the patient;
  - Prosthesis revision for hip instability;
- Given the considerations mentioned above, hip replacement is advised when the following conditions are met:
- The anatomical shapes and quality of the bone structures in patients must enable a prosthesis to be anchored;
  - The surgeon must inform patients of the risks of implanting a prosthesis and their agreement to the suggested operation must be obtained.

#### The following circumstances require special precautions:

- Patient is overweight;
- Excessive weight-bearing expected owing to profession or sports practices;
- Epilepsy or other factors that may cause repeated accidents presenting an increased risk of fractures;
- Osteoporosis and/or significant osteomalacia;
- Old infectious diseases with potential for relapse, with possibility of joint manifestations;
- Significant deformation of the joint to be replaced, indicating possible difficulties in anchoring the prosthesis;
- Intellectual incapacity of patients to understand and follow the surgeon's instructions;
- Alcoholism or drug addiction;
- Metabolic disorders or general drug treatments that may lead to progressive deterioration of the prosthesis bone support;
- Disorders in other joints (knee, ankle, etc.);
- Overall weak immune condition of the patient;
- Unbalanced diet.

#### 4 - CONTRAINDICATIONS

The following circumstances generally constitute an absolute contraindication to the implantation of a joint prosthesis:

- Allergy to any component of the prosthesis.
- Patients with a native cup that has a diameter less than 46 mm or greater than 62 mm.
- Evolving infection of the joint or relevant periparticular area.
- Loss of bone substances above or below the joint, making it impossible to anchor the prosthesis correctly;
- Severe neuroarthropathies; Charcot disease; muscle loss, neuromuscular lesion or vascular insufficiency of the affected limb, making the operation unjustifiable;
- Patients under 50 years old with joint lesions that could probably be treated successfully by a conservative treatment;
- Pregnancy.

#### 5 - IMPORTANT WARNINGS

##### 5.1 Preoperative warnings

- The operation must be scheduled with care. Depending on the results of the X-ray examination (templates) and preoperative planning, the surgeon will decide on the suitable implant.
- CERAVER® instruments and templates specially designed to prepare for the implantation of the joint components should be used. The description of the surgical procedure can never be exhaustive or fail to take into account all the risks and complications that must be accounted for.
- HAC cementless dual mobility cups must be used with a CERAVER® PE dual mobility insert of the same diameter.

##### 5.2 Postoperative warnings

- The joint prostheses must be handled with care. Scratches or unevenness due to impacts will lead to excessive wear and can give rise to complications.
- For the choice of prosthesis size: the templates, ancillary instruments and CERAVER® test implants will be the only ones used.
- Exposed prosthetic joint surfaces must not be in contact with any metal object, instruments for example.

##### 5.3 Anchoring procedure warnings

Anchoring is decisive in ensuring that the implant is reliably fitted. The following errors may give rise to mobilisation and complications:

- Excessive thinning of the cortical bone when preparing the bone bed.
- Selection inappropriate of the size of the prosthesis.
- Use of excessive force when fitting the prosthesis, leading to splintered fractures or a bone fracture line.

#### 6 - POSTOPERATIVE PRECAUTIONS

6.1 All patients implanted with a CERAVER® joint prosthesis should see an orthopaedic surgeon for postoperative check-up. This check-up is important to detect any signs of prosthetic wear or loosening before the appearance of functional manifestations. Early recognition of a complication enables effective prostheses to be taken in partnership with the available data on devices of similar composition, shape, and use. This evaluation has concluded there is conditional compatibility for MRI examinations up to 3 Tesla. However, as a precaution, it is recommended that MRI examinations are avoided in the 48 hours after implantation, and that the person responsible for the examination of the recent insertion of an implant is warned, if such an examination proves to be essential.

It should be noted that all devices presenting a high contrast with the surrounding biological environment may generate artefacts which must be taken into account for the perfect performance and interpretation of imaging examinations. For this purpose, a patient with an implant should be advised to inform the health professionals concerned (radiologists and radiology technicians) about the presence of an implant as soon as possible before these examinations.

#### 7 - COMPLICATIONS

The complications associated with cementless dual mobility cups are as follows:

##### Complications associated with the cups :

- Loosening;
- Luxation; Intra-prosthetic luxation;
- Periprosthetic fracture; cracks;
- Periprosthetic osteolysis; foreign body reactions;
- Traumatic fracture;
- Wear and tear, fracture of implant;
- Pinching (impingement);
- Infection of the prosthesis.

##### Complications associated with the procedure :

- Superficial and deep infections;
- Wound healing disorders;
- Claudication; Lower limb length inequality, Walking discomfort; Groin pain.
- Musculotendinous lesions;
- Cardiovascular complications Phlebitis; Deep vein thrombosis; Pulmonary embolism;
- Pulmonary complications;
- Dysplastic complications, paralytic ileus, intestinal necrosis;
- Neurological complication, nerve lesions, mental confusion;
- Other associated complications resulting from anaesthesia, use of ancillary equipment, nosocomial diseases

#### 8 - COMPATIBLE MEDICAL DEVICES

- Dedicated ancillary equipment.
- PE dual mobility liner with outer diameter appropriate to the inner diameter of the dual mobility cup

The CERAVER® PE dual mobility liners are compatible with CERAVER® 28 diameter stainless steel and alumina heads and all the CERAVER® stems (cemented Osteal and Ceratit) in accordance with the recommendations outlined in the respective instructions for each of these components.

#### 9 - STORING AND HANDLING PROSTHESES

- Joint prostheses must be stored in the original packaging.
- Joint prostheses must not be subjected to treatment or modification.
- Prostheses that are damaged or have undergone a treatment that is inappropriate or not authorized by the manufacturer must not be implanted.
- Do not use the device after the expiry date. Contact your CERAVER® representative and return the product to him/her.

#### 10 - PACKAGING AND STERILISATION

The various components are packed individually in double packaging and are supplied sterile. All CERAVER® hip implants are sterilised by gamma radiation at a minimum dose of 25 kGy. It is important to check the integrity of the sterility protection before use.

#### 11 - RESTERILISATION AND REUSE

Resterilisation and/or reuse of any of the components of a CERAVER® hip prosthesis is strictly prohibited due to the risk of exposing the patient concerned to contamination that could cause an infection.

#### 12 - DISPOSAL AFTER USE

If a device removed for revision surgery, dispose of the components according to the centre's medical waste disposal procedure. This does not apply if the follow-up procedure is due to a medical device vigilance case.

#### 13 - NOTE

The HAC cementless Dual Mobility Cup range was CE marked in 2004. Federal law in the USA and Canada requires that prosthetic components be purchased by a pharmacist. For further information concerning the use of these components, please contact the manufacturer, their overseas subsidiaries and distributors.

CERAVER - Les Laboratoires Ostéal Médical  
69 rue de la Belle Etoile  
95057 Roissy CDG Cedex - France  
Tél.: 33 (0)1 48 63 88 63 - Fax: 33 (0)1 48 63 88 99  
www.ceraver.com

Instructions last updated: February 2025

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