CERAFIXGENTA® BV

**STERILE R****STERILE EO****STERILE A**

CERAVÉR – Les Laboratoires Osteal Médical
69, rue de la Belle Etoile – 95957 ROISSY CDG cedex
Tél : +33 1 48 63 88 63 – Fax : +33 1 48 63 88 99
adv@ceraver.com – www.ceraver.com

● NOTICE D'INSTRUCTIONS

Français

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● INSTRUCTIONS FOR USE

English

Instructions relating to the following product:

Reference	Description
789	Low viscosity bone cement with gentamicin

Graphic symbols applicable to this product:

Symbol	Symbol title	Symbol	Symbol title
	Manufacturer		Do not reuse
	Expiry date		Sterile
	Batch code		Sterilised by irradiation
	Catalogue reference		Sterilised by ethylene oxide
	CE marking followed by the code for the notified body (0459, LNE/G-MED)		Sterilised by aseptic techniques
	Moisture sensitive		Do not resterilise
	Keep away from sunlight		Do not use if the packaging is damaged
	Do not keep at temperatures exceeding 25°C or 77°F		Warning
	Flammable		Consult instructions for use

Before using CERAFIXGENTA® BV cement, the surgeon should note the following:

1. PRODUCT DESCRIPTION

CERAFIXGENTA® BV cement is a cement based on radiopaque acrylic resins, specially designed for fixation of endoprostheses to be cemented to the bone. The functional characteristic of this product is that it is malleable for a short period, then hardens to retain a final shape, providing a stable, durable anchor for osteoarticular prostheses by evenly distributing the stress transferred from the cement to the bone when the patient moves. The cementation should, incidentally, be clearly visible on x-rays, thanks to the radiopacity conferred on the cement by the presence of zirconium dioxide.


This cement contains gentamicin to prevent colonisation of the prosthesis and surrounding tissues by germs sensitive to gentamicin.

This cement is obtained by mixing two pre-dosed, single-use, sterile components, one in liquid form, the other in solid form.


Liquid component:

With a volume of 21 ml contained in a brown glass ampoule, it has the following composition:

Methyl methacrylate	17.7 ml
N-butyl methacrylate	2.9 ml
N, N-dimethyl paratoluidine	0.4 ml
Hydroquinone	45 ppm

Methyl methacrylate


N-butyl methacrylate

Solid component:

With a volume of 48.0 g contained in a polyethylene sachet, it has the following composition:

PolyMethyl methacrylate	41.7 g
Zirconium dioxide	4.3 g
Activated benzoyl peroxide	0.5 g
Gentamicin in the form of sulphate equivalent to	0.80 g of gentamicin base

PolyMethyl methacrylate

When the liquid and solid components are mixed, the N, N-dimethyl paratoluidine breaks down the benzoyl peroxide to form free radicals, which allow polymerisation of the monomers to occur. This has the effect of binding the polymer microbeads together. Copolymerisation of the methyl methacrylate with a small quantity of n-butyl methacrylate improves the cement's impact resistance. The presence of zirconium dioxide enables cement to be detected on radiographs.

Gentamicin, an antibiotic from the aminoglycoside family, is heat-stable and water-soluble. Its bactericidal action derives from its ability to inhibit protein synthesis leading ultimately to cell lysis. The impact of this mode of action depends directly on the concentration of gentamicin. The minimum inhibitory concentration (MIC) is between 0.06 and 8 µg/mL.

Gentamicin is effective against many gram-negative bacteria (*Escherichia coli*, *Salmonella*, *Pasteurella* sp., *Pseudomonas aeruginosa*). Gram-positive bacteria are less sensitive, except for *Staphylococcus aureus* and *Staphylococcus epidermidis*, two important strains in the surgical field. Gentamicin is not active against Streptococcal species or anaerobic bacteria.

Combining it with cement produces bactericidal activity at the implant site. The *in vivo* results show bone levels significantly higher than the MIC after 18 months of implantation.

2. INDICATIONS

CERAFIXGENTA® BV is a radiopaque, low-viscosity, acrylic resin-based cement containing gentamicin, intended for the fixation of prostheses to be cemented to the bone in hip and knee replacements.

CERAFIXGENTA® BV cement is particularly recommended in re-operations for prosthetic infection.

3. CONTRAINDICATIONS

The following circumstances constitute a contraindication for the use of CERAFIXGENTA® BV cement:

- Allergies to the product components.
- Allergies to gentamicin or any other aminoglycoside.
- It is not recommended to use CERAFIXGENTA® BV cement for vertebroplasty or kyphoplasty procedures.
- Contraindication for fitting the prosthesis associated with cement

4. USE DURING PREGNANCY AND BREASTFEEDING

The safety of bone cement during pregnancy or breastfeeding has not been established. As a precaution, bone cement should not be used during the first three months of pregnancy and should be used during the rest of the pregnancy only in the event of a life-threatening disease.

5. CAUTIONS

- To ensure safe and effective use of CERAFIXGENTA® BV bone cements, the surgeon should have received specific training in their use and be familiar with the specifics of their handling, limits of use and application technique.
- As each component is pre-dosed, all of the solid component must be mixed with all of the liquid component.
- No product must be added or come into contact during preparation and application of this cement, to avoid any polymerisation problem or change in the physical and mechanical characteristics.
- Use of iodinated derivatives or rifamycin-class antibiotics can slow down or even inhibit the polymerisation reaction.
- As the liquid monomer is volatile and flammable, the operating theatre must be properly ventilated and all precautions relating to the use of flammable products must be taken.
- Any contact of the monomer with the skin, mucous membranes or latex items (surgical gloves) must be avoided.
- Prolonged inhalation of monomer vapours may cause drowsiness. Excessive exposure to concentrated vapours may also cause irritation of the respiratory tract and eyes.
- Concentrated monomer vapours may have adverse effects on soft contact lenses.

6. ADVERSE EVENTS

The following complications have been observed when using surgical cements based on methyl methacrylate:

- pulmonary embolism,
- trochanteric bursitis,
- vein thrombosis,
- heterotopic osteogenesis,
- haemorrhages and haematomas,
- cardiac arrhythmias and hypotension,
- stroke,
- reduced or low range of movement,
- fracture,
- rare allergic reactions to gentamicin,
- nephrotoxicity, ototoxicity: Gentamicin belongs to the aminoglycoside family. For this family, cases of renal failure and cochleovestibular damage have been reported in the event of high doses, prolonged treatment, prior renal impairment, haemodynamic disorders or combination with nephrotoxic and/or ototoxic products.

The surgeon is responsible for any complications that may arise from an erroneous indication and defective operating technique or incorrect asepsis; these cannot be attributed to the manufacturer or the supplier of the CERAFIXGENTA® products.

7. PREPARATION – INSTRUCTIONS FOR USE

- Preparation is carried out in the sterile area.
- The cement components should be at the ambient temperature of the operating theatre.
- Prepare the following equipment stored at operating theatre temperature:
 - a dry, sterile cup (polyethylene or ceramic or glass),
 - several sterile spatulas,
 - a cement syringe if the cement is to be injected using this system,
 - a timer.

7.1 MIXING PHASE

- Pour the whole of the solid component into the container and mix well with a spatula.
- Pour the whole of the liquid component over the powder. Start the timer.
- Work the liquid into the solid mass using the spatula.
- Mix thoroughly until a uniform texture is obtained, avoiding formation of air bubbles.

7.2 WAITING PHASE

SYRINGE

- Allow to stand, away from any heat source,
- Pour the liquid cement into the dedicated syringe for the application, avoiding formation of air bubbles.

MANUAL

- Allow to stand, away from any heat source

7.3 DELIVERY OF CEMENT PHASE

SYRINGE

- Inject the cement slowly, either from the top with the CERAVER® short nozzle syringe after placing a drain in the medullary canal, or from below with the CERAVER® long nozzle syringe.

MANUAL

- Take the paste without kneading when it no longer sticks to the gloves. It is then smooth.
- Apply the cement.

- Insert the implant
- Remove excess cement


7.4 HARDENING PHASE

- Keep the implant perfectly unmoving.

8. TIMELINE AND CLINICAL USE

Bone cements are heat sensitive. Any increase or decrease in temperature (ambient temperature and/or temperature of the cement components and the mixing equipment) will affect the handling characteristics and setting time of the cement. Variations in humidity are also associated with changing the handling characteristics and setting time.

Preparation chart

9. ASSOCIATED PRODUCTS

The table below shows the CERAVER® instruments associated with acrylic cements and their functions:

Surgical steps	Function	Reference	Designation
Preparation of surgical site	hold	615	Syringe holder
		616	Long-tip syringe rod
		614	Short-tip syringe stopper
	mix	799	Bowls and spatulas
Implant insertion	inject	621	Short-tip syringe
		619-PN	Long-tip syringe
		3549	Cement gun

Use of other instruments to apply CERAFIXGENTA® BV cement has not been validated by CeraVer®.

10. PRECAUTIONS FOR USE

- In order to prevent any contamination of the cement with fragments of glass, do not break the ampoule containing the liquid over the mixing apparatus.
- When used, the cement components and all the instruments with which they are in contact should be at operating theatre temperature.
- The following factors may affect the handling and setting time characteristics of the cement:
 - precise temperature of the components and instruments used for the preparation,
 - temperature and humidity level of the operating theatre,
 - preparation conditions,
 - shelf life.
- When using CERAFIXGENTA® BV for cementation of endoprostheses, it is essential to wash and vacuum clean all osteotomised surfaces that will be in contact with the cement immediately before application. After polymerisation, in order to avoid any contamination of articular surfaces with cement particles, it is imperative to clean the area surrounding the joint thoroughly.
- During and after introduction of the cement, the patient's blood pressure, pulse and breathing rate must be carefully monitored. A drop in blood pressure accompanied by a reduction in cardiac activity and acceleration of breathing rate may occur when the cement is introduced.
- Gentamicin nephrotoxicity and ototoxicity call for the following precautions for use:
 - In cases of renal failure or cochleovestibular damage, use gentamicin only if strictly necessary. Medical monitoring of kidney and hearing functions is required.
 - Avoid combining gentamicin with very active diuretics.
 - Do not administer cement with gentamicin at the same time as other potentially ototoxic and nephrotoxic substances.
- Gentamicin may potentiate the action of curare, myorelaxants and general anaesthetics. At the time of the procedure, the anaesthetist must be notified of the insertion of cement.
- As the cement polymerisation reaction is exothermic, the product reaches higher temperatures than the physiological level during the hardening phase.

11. STORAGE REQUIREMENT

- Cements must be stored in their original packaging.
- Cements must be stored in a dry, dark place and at a temperature lower than 25°C to prevent premature polymerisation of the liquid monomer.

12. STERILISATION

- The liquid component is sterilised by sterile filtration. The singly packaged ampoule is sterilised with ethylene oxide.
- The solid component is packaged in a polyethylene sachet, which is contained within a peelable outer sachet. The whole package is sterilised by gamma radiation (minimum dose of 25 kGy).
- Sterility is guaranteed only if the various packs of components are intact.

- The product will be sterilised only by the manufacturer. Regardless of the circumstances, no component may be resterilised.
- The duration of sterility is three years from the date of sterilisation, strictly on condition that the packaging remains intact. Do not use after the expiry date indicated on the packaging.

13. RESTERILISATION AND REUSE

Resterilisation and/or reuse of the components of a CERAFIXGENTA® BV cement unit (ampoule, sachet) is strictly prohibited because of the risk of exposing the patient concerned to contamination that could cause an infection. Reuse of the components of a CERAVÉR® bone cement would lead to a reduction in performance of the product concerned.

14. DISPOSAL AFTER USE

In the event of product returns, any disposal is managed by CERAVÉR®. If products are not returned, their disposal is the responsibility of the distributor or health care facility in accordance with local regulations.

15. DATE OF FIRST CE MARK

CERAFIXGENTA® BV low-viscosity orthopaedic cements were CE marked in 1996.

Federal law in the USA and CANADA requires that prosthetic components be purchased by a pharmacist.

For further information concerning the use of these medical devices, please contact the manufacturer, its overseas subsidiaries or distributors.

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