Réf. : IFU\_H26\_0814

## ● NOTICE D'INSTRUCTIONS

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### Français

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Femoral heads

### English

#### Instructions for use for the following implants:

- Alumina femoral heads (ref.: 2000, 2001, 2002, 2003, 2004, 2005)
- Stainless steel femoral heads (ref.: 2220, 2221, 2222, 2263, 2264, 2223, 2224, 2225, 2261, 2262, 2217, 2218, 2219)

Graphic symbols used for labelling these devices:

Symbol	Symbol title	Symbol	Symbol title
	Manufacturer		Do not re-use
	Expiration date		Device terminally sterilised by irradiation
	Batch code		N.B.: Not all the information required for use of the device is indicated on the label; it is necessary to consult the instruction manual.
	Catalogue reference of the device		

Prior to utilisation of CERAVER® hip joint prosthesis, surgeons should familiarise themselves with the following recommendations:

#### 1- GENERAL POINTS

- 1.1 The manufacturer and the supplier of the CERAVER® femoral head provide a warranty for the design and material quality of the product. A prosthesis which has already been used must not be implanted again. Instructions for assembling the different components must be respected, as defined in the commercial documentation. It is strictly forbidden to use CERAVER® hip replacement components together with implants of another origin.
- 1.2 The following factors are of prime importance for the success of an operation using a CERAVER® femoral head:
  - Careful selection of a suitable patient for the operation in question.
  - Exact surgical technique and proper use of special CERAVER® instruments.
  - Extremely strict asepsis.
- 1.3 The surgeon is responsible for any complications that may result from an erroneous indication, incorrect surgical technique or inadequate aseptic precautions. Neither the manufacturer nor the supplier of the CERAVER® implants can be blamed for these complications.
- 1.4 The surgeon must be completely familiar with all the aspects and problems of joint replacement surgery. The surgical procedures for the implantation of a CERAVER® prosthetic hip can be best learned at a clinic with many years' expertise in the use of these implants.
- 1.5 It should be noted that for patients receiving a hip joint prosthesis that their weight and activity level may affect the longevity of the prosthesis.

#### 2- DESCRIPTION OF IMPLANTS

CERAVER® femoral heads are designed and manufactured in keeping with the latest developments in technology. Their mechanical and functional characteristics can only be guaranteed if the surgical technique and the present instructions are strictly observed.

- The femoral heads are made of FeCrNiMnMo stainless steel with a high nitrogen content (ISO 5832-9 standard, ASTM F 1586 standard), or Al<sub>2</sub>O<sub>3</sub> ceramic (ISO 6474-1 standard).

#### 3- INDICATIONS

- A joint prosthesis should be considered only after all other surgical methods of treatment have been carefully weighed up against each other and none have been judged to be more appropriate.
- The most successful implants are those which are implanted in a sound natural joint. On the other hand, an artificial joint can be a highly beneficial substitute for a deformed or diseased joint, because it eliminates pain and is conducive to the restoration of good mobility and weight-bearing capacity.
- Every artificial joint is always subject to wear, which still remains a major problem that cannot be avoided. An initially stable prosthesis may become loose in the course of time. Wear and loosening of a prosthesis may make revision surgery necessary.
- Infection of the artificial joint may have a negative effect on most patients. Usually it necessitates removal of the prosthetic components and may under certain circumstances lead to incapacitating residual effects.

#### Planned applications:

Hip arthroplasty consists of fully or partially replacing a hip joint suffering from various affections that cannot be treated through other therapies or methods. It aims to restore all or part of the patient's mobility and stability while minimising pain or even eliminating it. The replacement of the hip joint may be indicated in the following cases:

- Primary coxarthrosis;
- Osteonecrosis of the femoral head;
- Proximal femoral fracture;
- Post-traumatic arthritis;
- Sequelae following injury;
- Non-infectious arthritis;
- Legg-Calve-Perthes disease;
- Epiphysiolysis;
- Hip dysplasia;
- Sequelae of hip infections;
- Dislocated hip;
- Gaucher's disease;
- Arthrodesis;
- Total hip replacement revision.

In general, the patient must be suffering from severe arthropathy that does not respond to any conservative therapy or other alternative treatment. Based on the above statements, a joint prosthesis is indicated in cases where the following conditions are fulfilled:

- The selected patient's joint is anatomically and structurally suited for the implantation of a prosthesis.
- Patients must have been informed by the surgeon of the risks associated with the implantation of a joint prosthesis and have declared their consent to the operation.

Comment: in the event of revision surgery following fracture of the alumina head, given the presence of numerous particles in the tissue surrounding the implant, it is preferable that the cone of the cone present -- to re-implant an alumina head in order to avoid any risk of abrasion of the metal head.

The following circumstances call for particular caution:

- Overweight patients;
- Excessive loading to be expected as a result of work and sport;
- Epilepsy or other factors favouring repeated accidents with increased risk of fracture;
- Severe osteoporosis and/or osteomalacia;
- Past history of sites of potentially active infectious diseases with potential joint manifestations;
- Severe deformation of the affected joint, which may render fixation of the implant more difficult;
- Weakening of the supporting structures due to tumours;
- Mental inability of patients to understand and follow the surgeon's instructions;
- Drug or alcohol addiction;
- Metabolic disorders or general medical treatments, leading to a progressive deterioration of bone tissue;
- Other joint affections (knee, ankle, etc.);
- Patient's overall immune status weakened;
- Unbalanced diet.

#### 4- CONTRAINDICATIONS

The following conditions are generally accepted as contraindications for the implantation of a joint prosthesis:

- Allergy to the implant materials of which the prosthetic components are made;
- Active infection of the joint or surrounding area;
- Severe neuroarthropathy;
- Loss of bone structure proximally and distally of the joint, which precludes proper anchorage of the prosthesis;
- Significant muscle loss making joint prosthesis stability impossible, neuromuscular damage or vascular insufficiency of the limb affected, making the surgery unjustifiable;
- Charcot disease;
- Patients aged less than 50 years old and suffering from joint damage that could probably be treated successfully with a conservative treatment;
- Pregnancy.

#### 5- IMPORTANT WARNINGS

##### 5.1 Preoperative warnings

- The operation should be planned carefully. On the basis of the radiological templates and pre-operative planning, the surgeon will choose the appropriate implant.
- The special CERAVER® templates and ancillary instrumentation should be used to prepare the fitting and setting of the articular components. A description of the operative technique can never be exhaustive or contain all the risks and complications that should be considered. Technical product literature and the literature on the use of the implants, as well as videos, are available on request from CERAVER® France (or from your local CERAVER® representative).
- The male cone of the femoral stem has to present the same diameter than the female cone of the femoral head. The size of the cone is indicated on the box and on the implant.
- Femoral heads made of Alumina ceramic must be assembled with stems of the same origin that have never been used.
- Stainless steel femoral heads must not be used with the alumina ceramic liner of the CERAFIT cup.
- In the event of fractures to the alumina head, the surgeon must decide whether to replace the femoral head. The most appropriate solution, recommended by CERAVER®, is the removal of the femoral stem.

##### 5.2 Intraoperative warnings

- Joint prostheses have to be handled with care. Even slight scratches or impact dents on the articular surface may contribute to excessive wear and initiate undesirable complications.
- Only CERAVER® ancillary instrumentation should be used for selecting and fitting CERAVER® prosthetic components.
- Contact articular surfaces must not be marked or come into contact with metallic or other hard objects (e.g. instruments, table tops, etc.).
- Ceramic femoral heads should only be implanted on a perfectly clean and dry cone.
- In the event of fractures to the alumina head, the surgeon must decide whether to replace the femoral head. The most appropriate solution, recommended by CERAVER®, is the removal of the femoral stem.

#### 6- POSTOPERATIVE PRECAUTIONS

- 6.1 All patients with an articular CERAVER® prosthesis must be postoperatively monitored by an orthopaedic surgeon. Follow-up examinations are essential to detect signs of wear or prosthesis loosening before any resulting detrimental effects. Early detection of an impending complication allows the surgeon to initiate timely and effective countermeasures. Above all, a re-operation performed at the right time has a much better chance of success than a later surgical intervention, when advanced loosening has already resulted in serious resorption of the bone stock.
- 6.2 Patients must be instructed to inform their surgeon of the slightest changes in the joint operated on.
- 6.3 An annual clinical check-ups are not possible, the surgeon should be sent a control radiograph of the joint in question. This will help in the early detection of any impending complication allow the surgeon to initiate timely and effective countermeasures. Above all, a re-operation performed at the right time has a much better chance of success than a later surgical intervention, when advanced loosening has already resulted in serious resorption of the bone stock.
- 6.4 The innocuousness and compatibility, as well as overheating and migration of CERAVER® femoral heads have not been evaluated using an MRI.

#### 7- COMPLICATIONS

- The following complications have occasionally been observed and consequently require special attention:
  - Implant loosening resulting from changed load transmission conditions, or attrition and breakage of the cement bed and/or tissue reactions to the presence of the prosthesis. Loosening is frequently due to one or several risk factors listed above;
  - Fracture of the alumina head;
  - Early and late infections;
  - Dislocation, subluxation, inadequate mobility, undesirable shortening or lengthening of the limb operated on following non-optimal implantation of the prosthesis;
  - Excessive medialisation of the joint causing pain in other joints or in the back;
  - Bone fractures following unilateral excess loading or bone fatigue;
  - Transient or permanent nerve lesions due to compression or haematoma;
  - Cardiovascular complications or diseases, including phlebotrombosis, pulmonary embolism and Charcot disease. Damage to the blood vessels during the surgical procedure;
  - Wound haematoma and slow wound healing;
  - Tissue reactions as a result of allergy to the implanted material, in particular metal, or due to the accumulation of particles;
  - Temporary or permanent peripheral neuropathy;
  - Lung disorders, including in particular pneumonia and atelectasis;
  - Urological complications (retention, urinary tract infection);
  - Other complications arising from anaesthesia, the use of ancillary equipment, nosocomial infections;
  - Ectopic ossification;
  - Pain.

#### 8- COMPATIBLE MEDICAL DEVICES

- Stainless steel femoral heads are compatible with standard and anti-dislocation polyethylene liner with the same joint diameter.
- Alumina and stainless steel femoral heads are compatible with the following devices:
  - Dedicated ancillary instrumentation;
  - PE cemented cups with the same diameter;
  - Dual mobility PE inserts (only with 28 mm diameter femoral heads);
  - Bipolar heads;
  - CERAFIT cementless stems;
  - OSTÉAL cementless stems.

Comment: articular diameter means nominal diameter of bearing surface.

#### 9- STORAGE AND HANDLING OF PROTHÈSES

- Prosthetic components should be stored in their unopened original containers.
- Joint prostheses must not be subjected to treatment or modification.
- Prostheses that are damaged or have undergone a treatment that is inappropriate or not authorised by the manufacturer must not be implanted.
- Do not use the device after the expiration date. Contact your CERAVER® representative to return the product.

#### 10- PACKAGING AND STERILISATION

The various components are packed individually in double packaging and are supplied sterile. The sterilisation and/or reuse of prosthetic components is strictly prohibited under penalty of exposing the concerned patient to contamination that can cause infection.

#### 12- DISPOSAL AFTER USE

If a device is removed for revision surgery, please dispose of the components according to the centre's medical waste disposal procedure. This does not apply if the revision procedure is due to a Medical Device Vigilance case.

#### 13- NOTE

- The stainless steel femoral heads were granted EC marking in 1998.
- The alumina femoral heads were granted EC marking in 1998.
- Federal law in the USA and CANADA requires that prosthetic components be purchased by a pharmacist.
- For further information concerning the use of these components, please contact the manufacturer or its foreign subsidiaries and distributors.

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Date of last revision of this document August 2014

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