● NOTICE D'INSTRUCTIONS

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● INSTRUCTIONS FOR USE

Hip, Knee and Biomaterials Surgical instruments Recommendations for reprocessing reusable

1. PURPOSE

These instructions are recommended for the care, cleaning, disinfection, maintenance and sterilisation of Ceraver® reusable hand-held orthopaedic surgery instruments. This document has been created to assist health care personnel to achieve safe handling, as well as effective maintenance and reprocessing of Ceraver® reusable instruments.

2. SCOPE

This instruction manual provides information on the care, decontamination, cleaning, disinfection, maintenance and sterilisation of reusable surgical instruments manufactured and/or distributed by Ceraver®.

3. GLOSSARY

- **Contaminated:** state of an item having actually or potentially been in contact with micro-organisms.
- **Decontamination:** use of physical or chemical means to remove, deactivate or destroy blood-borne pathogens on a surface or item to the point where they can no longer transmit infectious particles and the surface or items are safe for handling or disposal.
- **Disinfection:** process used to reduce the number of viable micro-organisms on a product to a level previously specified as appropriate for handling or use.
- **Maintenance:** visual and/or functional inspection
- **Packaging:**
 - **Individual:** a standard packaging material may be used. Ensure that the bag is large enough to contain the instruments without pulling on the closure.
 - **Batch:** the instruments may be placed in instrument trays provided for the purpose or in general-purpose sterilisation trays. Ensure that the cutting edges are protected.
- **Washer/Disinfectant:** machine designed to clean and disinfect medical instruments and other items used in medical, dental, pharmaceutical and veterinary practice settings.
- **Cleaning:** removal of any contamination on an item, required before any other processing
- **Manual cleaning:** cleaning without using an automatic washer or washer/disinfectant

Note: cleaning and disinfection are often carried out in a single operation.

- **Sterile:** free from all viable micro-organisms.
- **Sterilisation:** validated process used to remove all forms of viable micro-organisms from an instrument.

Note: In a sterilisation process, the nature of microbiological death is described by an exponential function. Therefore, the presence of micro-organisms on an item can be expressed in terms of probability. Although this probability can be reduced to a very low number, it can never be reduced to zero. This probability can only be assured for validated processes.

- **Processing/Reprocessing:** Process consisting of decontamination, cleaning, infection, packaging and sterilisation, required to prepare a new or used surgical instrument prior to its use.

4. ACRONYMS

- BI = Biological Indicator
- PPE = Personal Protective Equipment
- TSE = Transmissible Spongiform Encephalopathy
- CJD = Creutzfeldt-Jakob Disease
- SAL = Sterility Assurance Level
- SSD = Sterile Supply Department
- OT = Operating Theatre
- AORN = Association of periOperative Registered Nurses

5. CONSIDERATIONS

This instruction manual refers to all Ceraver® reusable surgical instruments and must be read carefully. These instruments are designed specifically for the placement of Ceraver® orthopaedic implants (cf. table of implants §10).

Orthopaedic surgery requires heavy instruments comprising several components, articulated or rotating parts and removable handles. Instruments are often supplied in sets and divided into cases in which they are arranged by size or in the order required for a specific surgical procedure.

During musculoskeletal surgery, instruments are contaminated with blood, bone and tissue debris. They can also be contaminated with body fluids containing the hepatitis virus, HIV or other aetiological agents and pathogens. All health care personnel should be familiar with the universal precautions necessary to prevent injuries caused by handling instruments with cutting edges and/or sharp points during and after surgical procedures or during reprocessing.

These Ceraver® devices must only be used by medical health care personnel, such as surgeons, IBODEs (state-certified operating theatre nurses), qualified operating theatre nurses and reconditioning personnel (disinfection, washing and sterilization)/ They are trained in the handling and in mastering their use and the operating technique for the implant in question.

New and used instruments must be thoroughly processed according to these instructions prior to use.

It is the responsibility of healthcare facilities to clean, disinfect, package and sterilise all instrument sets before returning them to Ceraver® or to recondition them for a new use if the material remains in storage.

The user who receives Ceraver® equipment must also inspect it to verify that the instruments have been properly cleaned and decontaminated before repeating the reprocessing procedures to prepare the loan set for subsequent reuse.

Core orthopaedic instrument sets must be complete and in good condition to be used correctly. Optional instruments may be available on request from your Ceraver® representative.

To maintain the instruments correctly, it is important to consider the following information and processing instructions:

- Warnings and precautions
- Instrument set completeness and functionality
- Reprocessing limitations and/or restrictions
- Preparation for processing at the point of use
- Preparation for cleaning (including assembly/disassembly, if necessary)
- Cleaning, disinfection and drying
- Maintenance, inspection, testing and lubrication
- Sterile packaging
- Sterilisation
- Storage

6. PROCESSING INSTRUCTIONS

These surgical instruments are supplied non-sterile and must be cleaned, disinfected and sterilised prior to use. After use, these instruments must be decontaminated, cleaned, put away, packaged and sterilised. The following information describes the correct steps, based on Ceraver®'s experience and the recommendations of organisations (see references), for the processing of Ceraver® surgical instruments to ensure they have a safe, durable and efficient service life.

The instructions below describe the necessary processing steps that new and used instruments must undergo to attain sterility.

Warnings and precautions:

- **Universal precautions must be observed** by hospital personnel who work with contaminated or potentially contaminated medical instruments. Instruments with sharp points or cutting edges should be handled with extreme caution.
- **Personal Protective Equipment (PPE) should be worn** when handling or working with contaminated or potentially contaminated equipment, instruments and materials. PPE includes lab coat, mask, protective goggles or face shield, gloves and shoe covers.
- **Do not use metal brushes or scouring pads** during manual cleaning procedures. These materials may damage the surface and finish of the instruments. Use soft-bristled, nylon brushes and swabs.
- Cleaning agents with low foaming surfactants should be used during manual cleaning procedures to ensure that instruments are visible in the cleaning solution. Manual scrubbing with brushes should always be performed with the instrument immersed in the cleaning solution to prevent generation of aerosols and splashing, which may spread contaminants. Cleaning agents must be gently and thoroughly rinsed from instrument surfaces to prevent accumulation of detergent residues.
- **Do not place heavy instruments on top of delicate instruments.**
- **Do not allow contaminated instruments to dry before reprocessing.** All subsequent cleaning and sterilisation steps are facilitated by not allowing blood, body fluids, bone and tissue debris, saline solution or disinfectant to dry on used instruments.
- Saline solution and cleaning/disinfection agents containing aldehyde, mercury, active chlorine, chloride, bromine, bromide, iodine or iodide are corrosive and **must not** be used. Instruments **must not** be placed or immersed in **Ringer's solution**.
- Only instruments manufactured and/or distributed by Ceraver® should be included in Ceraver® instrument cases. These validated reprocessing instructions are not applicable to Ceraver® cases that include instruments that are not manufactured and/or distributed by Ceraver®.
- Descaling agents that include morpholine should not be used in steam sterilisers. These agents leave residues that can damage polymer instruments.
- It is important to take special care when manually cleaning instruments for boring, reaming or drilling so that all bone residue is properly removed.
- If an instrument breaks in situ, a fragment can be located using imaging techniques tailored to the material.

Inspection on receipt:

- Upon receipt, instrument sets should be inspected to check for completeness. In addition, illustrated bills of materials are provided with the cases containing stored instruments to check the presence of instruments on arrival and departure from the facility.
- Orthopaedic surgical procedures follow a specific order in which the instruments are used. In addition, many instruments have dimensional features which determine bone resections and implant sizes, and measure intramedullary canal sizes, depth of drill holes, angles of tubes/ plates, acetabular cup placements, etc. It is therefore very important that all required sizes of a specific instrument series are available (specific instruments are routinely omitted from instrument sets due to infrequent use, unless requested by the user). Contact your Ceraver® representative if requested instruments have been omitted but are required for the operation.

Limitations and restrictions:

Neutral pH enzymatic cleaning agents are recommended for cleaning Ceraver® reusable instruments. Alkaline agents with a pH of 12 or less can be used to clean stainless steel and some polymer instruments in countries where required by law or local regulations or where there is a risk of prion disease, such as Transmissible Spongiform Encephalopathy (TSE) and Creutzfeldt-Jakob Disease (CJD). **It is extremely important to neutralise and rinse alkaline cleaning agents from the instruments thoroughly and completely.**

Note: Drill bits, reamers, rasps and other cutting instruments should be carefully inspected after processing with alkaline detergents to ensure that cutting edges are fit for use.

Note: It is important to select enzymatic solutions intended for breakdown of blood, tissues and body fluids. Some enzymatic solutions are designed specifically for breakdown of faecal matter or other organic contaminants and may not be suitable for use with orthopaedic instruments.

- Repeated processing in accordance with the instructions in this manual has no degrading effect on the expected performance on Ceraver® reusable instruments unless otherwise specified. End of life for stainless steel or metal surgical instruments is determined by wear and deterioration due to the intended surgical use and not to reprocessing.
- Ceraver® does not state the maximum number of reuses permitted for ancillary equipment. The service life of reusable surgical instruments depends on various factors, including the kind of handling endured and the type of processing between uses. The best way to determine when a device should cease to be reused is to carry out a visual inspection and to check that is functioning properly before each use.
- At the end of the operation, the surgeon should take the utmost care not to leave instruments or instrument debris in situ, for the safety of the patient.
- Automatic cleaning using a washer/disinfector alone may not be effective enough for orthopaedic instruments with lumens, slots, cannulas, blind holes, contacting surfaces and other complex features. A manual cleaning process is recommended. This manual cleaning procedure can be combined with automatic cleaning, but automatic cleaning alone is not sufficient and is not recommended.**

Note: only manual cleaning has been validated by Ceraver®.

- Where applicable, multi-component instruments should be disassembled for cleaning. Disassembly, where necessary, is generally self-evident. Take care not to lose any parts. If a part is lost, notify your Ceraver® representative as soon as possible.
- Instruments must be removed from the metal or polymer cases for manual and/or manual/automatic cleaning procedures. Do not clean instruments while they are in the polymer or metal cases. Instrument cases and lids must be cleaned separately from the instruments.
- The polymers used in Ceraver® instrument sets can be sterilised using steam/moist heat. If polymer surfaces become «chalky», show excessive surface wear (e.g., cracking or delamination) or if the polymer instruments show visual or functional abnormality, they must be replaced. Notify the Ceraver® representative if polymer instruments need to be replaced.
- Instruments must not come into contact with fluorinated or chlorinated products, or with oily detergents. Instruments containing components made of synthetic materials (plastics), must on no account come into contact with strong acids (pH <4), bases, organic or ammonium solvents, or with oxidising chemicals or any other product that may damage the material.
- Steam/moist heat is the recommended sterilisation method for Ceraver® instruments
- Ethylene oxide (EO) sterilisation, plasma sterilisation and dry heat sterilisation methods are not recommended for sterilisation of Ceraver® reusable instruments.
- Instruments with removable polymer sleeves must be disassembled for sterilisation.
- Avoid using acidic anti-corrosion agents.
- Use of hard water should be avoided. Softened tap water may be used for initial rinsing. Purified water should be used for final rinsing to eliminate mineral deposits on the instruments. One or more of the following processes may be used to purify water: ultrafiltration (UF), reverse osmosis (RO), deionisation or equivalent.

Point of use – preparation for processing:

- Remove excess body fluids and tissue from the instruments with a lint-free, disposable cloth. Place instruments in a basin of distilled water or in a case covered with a damp cloth. Do not allow saline, blood, body fluids, bone and tissue fragments or other organic debris to dry on instruments prior to cleaning

Note: soaking in proteolytic enzyme solutions facilitates cleaning, especially for complex instruments with hard-to-reach areas (e.g., tubular designs, etc.). These enzymatic solutions dissolve proteinaceous matter and prevent blood and protein-based materials from drying on the instruments. The manufacturer's instructions for preparation and use of these solutions must be followed explicitly.

- The instruments should be cleaned within 30 minutes of use to minimise the possibility of drying prior to cleaning.
- Used instruments must be transported in closed or covered containers to prevent any unnecessary contamination risk.

Preparation before cleaning:

- Symbols or specific instructions engraved on the instruments or on the instrument cases should be strictly followed.
- Where applicable, multi-component instruments should be disassembled for appropriate cleaning. Care should be taken to avoid losing small screws and components. If a part is lost, notify your Ceraver® representative as soon as possible.
- Published instructions for use and surgical techniques and/or procedures may serve as an additional source of information to illustrate the assembly/disassembly instructions for specific Ceraver® instruments.

Preparation of cleaning agents:

- Ceraver® recommends using neutral pH enzymatic cleaning agents with low foaming surfactants. Alkaline agents with a pH of 12 or less can be used in countries where the laws or local regulations require it. Use of alkaline agents should be followed by a neutralising agent and thorough rinsing.
- All cleaning agents should be prepared in accordance with the dilution and temperature recommendations of the manufacturer. Softened tap water may be used to prepare cleaning agents. Use of the recommended temperatures is important for optimal performance of cleaning agents.
- Dry powder cleaning agents should be completely dissolved prior to use to prevent oxidation or corrosion of the instruments.
- Fresh cleaning solutions should be prepared when existing solutions become visibly contaminated (bloody and/or cloudy).

CLEANING PRODUCTS	
Water	Use cold demineralised water or reverse osmosis filtered water. Temperatures higher than 60°C (140°F) should be avoided, as they have a coagulating effect on proteins, making them difficult to remove.
Detergent	Prepare a detergent (e.g., LIQUI-NOX®, Alconox, Inc. pH 8.5) according to the manufacturer's instructions.
Enzymatic cleaner	Prepare an enzymatic cleaner (e.g., ENDOZIME®, Ruhof Corporation pH 6.0 to 7.5) according to the manufacturer's instructions.
Accessories for manual cleaning	Soft nylon bristle brushes and/or bottle brush, syringes, gloves, disposable absorbent cloths (e.g., KIMWIPE®, Kimtech Science)
Ultrasonic cleaner	Ultrasonic cleaners must be checked regularly to ensure that they are functioning correctly

Instructions for manual disinfection/cleaning:

CLEANING/DISINFECTION		
Manual cleaning	Step 1	Immerse the instruments completely in an enzymatic solution and allow to soak for 20 minutes at room temperature (+15°C / +25°C).
	Step 2	Scrub using a soft nylon bristle brush for four minutes (or five minutes for large items) in the enzymatic solution at room temperature (+15°C / +25°C). Add a rinse for items with holes, lumens and hard-to-reach areas.
	Step 3	Remove instruments from the enzymatic solution and rinse with tap water for at least three minutes. Rinse lumens, holes and all other difficult-to-reach areas carefully and thoroughly.
	Step 4	Put the prepared cleaning agents inside a sonication unit. Immerse the instrument completely in the cleaning solution and sonicate for 10 minutes at 45-50 kHz at room temperature (+15°C / +25°C).
	Step 5	Scrub using a soft nylon bristle brush for four minutes (or five minutes for large items) in the enzymatic solution at room temperature (+15°C / +25°C). Add a rinse for items with holes, lumens and hard-to-reach areas.
	Step 6	Rinse the instrument twice in demineralised water for at least three minutes each time at room temperature (+15°C/+25°C). Rinse lumens, holes and all other difficult-to-reach areas carefully and thoroughly.
	Step 7	Remove excess moisture from the instrument with a clean, absorbent, lint-free cloth.

These manual disinfection and cleaning instructions can be followed by automatic cleaning:

CLEANING STEPS	
Manual and automatic cleaning	1. Repeat steps 1 to 7 of the manual cleaning process 2. An automatic cleaning procedure can involve a washer-steriliser, a washer-disinfector, an ultrasonic washer or other types of machines for washing and decontaminating instruments. There are many types of automatic washer systems, each with their own instructions that must be followed. 3. Washer-disinfector systems are not recommended as the sole cleaning method for surgical instruments. Orthopaedic instruments should be cleaned following the manual cleaning procedure or the combined manual/automatic cleaning procedure outlined in this manual, unless otherwise specified. 4. An automatic washer/disinfector can be used following the above-mentioned manual cleaning procedure, but it is not compulsory. 5. Simple instruments without multiple components, lumens/slits/cannulas, blind holes, contacting surfaces, connectors and internal mechanisms or other complex features may be cleaned and disinfected using a standard washer/disinfector cycle for surgical instruments. Instruments should be thoroughly inspected prior to sterilisation to ensure effective cleaning.

Inspection, maintenance, testing and lubrication

1.

Carefully inspect each instrument to ensure that all visible contamination has been removed. If contamination is found, repeat the cleaning/disinfection process.
2.

Visually inspect for completeness, absence of damage and/or excessive wear.

Note: If there is damage or wear that could compromise the function of the instrument, contact the Ceraver® representative for a replacement.

3.

Check that moving parts are functioning properly. (e.g., joints, hinges, sliding parts, etc.).
4.

Hinged, rotating or articulated instruments should be lubricated with a water soluble product (e.g., instrument milk or an equivalent lubricant) intended for surgical instruments to be sterilised.

Note: Do not use mineral oil or silicone lubricants, as they:
a. coat micro-organisms;
b. prevent direct contact of steam with surfaces;
c. are difficult to remove.

Note: These instructions do not apply to electric or pneumatic instruments. These instruments have different requirements and should be lubricated according to the manufacturer's instructions.

- 5. Check elongated instruments (particularly rotating instruments) for deformation.
- 6. For instruments that need to be assembled, check that they fit together properly: Mixing instruments or components from different suppliers is not recommended, for mechanical and design compatibility reasons. CERAVÉR accepts no responsibility for incompatibility when instruments or components from different sources are used. This recommendation does not apply to instruments with universal connections (AO-type connections).

Packaging:

- **Packaging instrument sets individually:**
- Cleaned and disinfected instruments should be packaged in their disassembled state.
- When packaging individual instruments, use single use medical quality pouches suitable for steam sterilisation and of appropriate size, and double wrap the instruments.
- **Instrument cases and trays with defined, preconfigured layouts should be used under the following conditions:**

- Cleaned and disinfected instruments should be packaged in their disassembled state.
- When putting instruments into sterilisation cases with lids, use a medical-grade paper suitable for steam sterilisation and wrap the cases using the double-wrapping method. Cases with lids can also be placed in an approved sterilisation container with a sealed lid for sterilisation.
- Wrapping and packaging must comply with standard EN ISO 11607, and must also be suitable for steam sterilisation (temperature resistance up to at least 138°C (280°F), with sufficient steam permeability).
- Areas designated for specific instruments must contain only instruments specifically intended for these areas.
- Optional Ceraver® instruments should not be added to a preconfigured instrument case unless a dedicated universal space or compartment has been included in the design and the guidelines described below for cases without defined layouts or universal spaces can be applied.
- Only instruments manufactured and/or distributed by Ceraver® should be included in Ceraver® instrument cases. These validated reprocessing instructions are not applicable to Ceraver® cases that include instruments that are not manufactured and/or distributed by Ceraver®.

- **Universal instrument cases and trays without defined, preconfigured layouts or containing undefined universal spaces or compartments should only be used under the following conditions:**

- If possible, all instruments must be disassembled before being placed in the case.
- Instructions for assembly and disassembly are included with the corresponding bill of materials and sent to the operating theatre for assembly/disassembly of disassembled sterilised instruments.
- All the instruments should be arranged to ensure steam penetration to all instrument surfaces. Instruments should not be stacked or placed in close contact.
- The user must check that the instrument case is not tipped or the contents shifted once the devices are arranged in the case.
- Only instruments manufactured and/or distributed by Ceraver® should be included in Ceraver® instrument cases. Ceraver® validated reprocessing instructions are not applicable to Ceraver® cases that include instruments that are not manufactured and/or distributed by Ceraver

Sterilisation instructions:

- Refer to the table below for the recommended minimum sterilisation parameters that have been validated by Ceraver® to offer a Sterility Assurance Level (SAL) of 10-6.

Cycle type	Minimum temperature	Minimum exposure time	Minimum dry time
Prevacuum	134°C	18 min	30 min

Note: the steriliser manufacturer's instructions for operation and load configuration should be followed explicitly.

- The facility is responsible for the internal procedures for the reassembly, inspection and packaging of instruments once they have been thoroughly cleaned in a manner that will ensure steam sterilant penetration and adequate drying. Provisions for protection of any sharp or potentially dangerous edges of the instruments should also be recommended by the care facility. In all cases, because of the differences between sterilisation devices and the microbial load on the appliances when they are used in clinic, it is generally the responsibility of each user to check and verify that the sterilisation cycle used will give an appropriate sterility assurance level for the instruments.
- Steam/moist heat sterilisation is the recommended method for Ceraver® orthopaedic instrument sets.
- Instrument sets should be properly prepared and packaged in Ceraver trays designed to allow steam to penetrate and make direct contact with all surfaces.
- Ethylene oxide or plasma sterilisation methods should not be used.
- Gravity displacement sterilisation cycles are not recommended because cycle times are too long.

Storage instructions:

- Containers containing instruments that are sterile and/or wrapped in double bags should be stored in a designated limited access area that is well ventilated and protected from dust, moisture, insects and extreme temperatures/humidity.
- Sterile instrument packages should be carefully examined prior to opening to ensure that package integrity has not been compromised.

Note: maintenance of sterile package integrity is generally event related. If a sterile wrap is torn or perforated (in the case of «draping»), shows evidence of tampering or has been exposed to moisture, the instrument set must be cleaned, repackaged and sterilised.

Note: if there is any evidence that the lid seal or filters on a sterilisation container have been opened or handled, the sterile filters must be replaced and the instrument set resterilised.

7. FACILITIES' RESPONSIBILITIES FOR SETS OF CERAVÉR® INSTRUMENTS ON LOAN OR IN DEPOSIT

- Orthopaedic surgical instruments generally have a long service life; however, mishandling or inadequate protection can rapidly reduce it. Instruments which no longer perform properly because of long use, mishandling or improper maintenance should be returned to Ceraver®. Notify your Ceraver® representative of any problem with the instruments.

- Loan sets should undergo all the steps of decontamination, cleaning, disinfection, inspection and final sterilisation before being returned to Ceraver®. The shuttle sheet documenting the decontamination and sterilization procedure should be provided with instruments being returned to Ceraver®.
- Missing or damaged instruments from loan sets should be brought to the attention of the operating theatre referent to forward to the Ceraver sales administration and your Ceraver® representative to ensure the traceability of the event, the any costs incurred and the information to the Ceraver reception service so that the restoration of the complete instrument set is working order for the next user.
- The instructions provided in this manual have been validated by Ceraver® and are fully tailored to the preparation of orthopaedic devices for use. It is the responsibility of the facility to ensure that reprocessing is performed using the appropriate equipment and materials, and that the personnel tasked with the reprocessing have been correctly trained in order to achieve the desired result. The equipment and processes should be validated and monitored regularly. Any deviation from these instructions by the processing supervisor should be properly assessed for effectiveness to avoid any potential adverse consequences.

8. CLAIMS AND REQUESTS FOR INFORMATION

Ceraver® makes every effort to design and supply devices that are safe and effective, in accordance with its Quality policy. However, in the event of a serious incident linked to a Ceraver® device, the health professional (user, prescriber, etc.) or the patient should contact Ceraver® immediately either by the Ceraver® representative of the sector or the Quality department or the Sales administration.

Claims or any other reasons for dissatisfaction concerning Ceraver® devices should be stated, specifying the reference, batch number and a full description of the incident.

For more information, please contact your Ceraver® representative or customer service directly using the contact details below:

CERAVÉR - Les Laboratoires Osteal Médical
69 rue de la Belle Etoile
95957 Roissy CDG Cedex - France
Tél. : 33 (0)1 48 63 88 63 - Fax : 33 (0)1 48 63 88 99
www.ceraver.com - ADV@ceraver.com
Date of last review: September 2023

9. REFERENCES

1. ISO 17664, Sterilisation of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices
2. ISO 17665-1, Sterilisation of Health Care Products – Moist heat

Compliance with standards ISO 17664 and ISO 17665-1 is noted in the sterility and cleaning validation reports.

10. NOTE

The CE mark dates for ancillary equipment are set out below.

Ancillary equipment for	Dates of first CE mark
Shoulder	
ULYS total anatomic shoulder prosthesis	2006
Knee	
Tricompartmental knee prosthesis	
HERMES PS REVISION range	1997
HIFIT PS range	2012
Cemented unicompartmental knee prosthesis	
HERMES FP	1998
HERMES UNI	1990
Additional components	
Tibial extension	1998
Patellar button	1996
Hip	
CERAFIT – HAC Femoral stems	
CERAFIT R-MIS HAC Ilateralised stem	2011
CERAFIT R-MIS R stem	2004
CERAFIT Multicone HAC stem	1997
OSTEAL – Cemented femoral stems	
OSTEAL standard stem	1977
OSTEAL revision stem	2004
Bipolar head	1996
Alumina femoral head	1998
Stainless steel femoral head	1998
Alumina insert	2004
Polyethylene dual mobility insert	2004
PE insert	2004
PE cemented cup	1995
CERAFIT cup screw	1995
HAC dual mobility cup	2004
Cemented dual mobility cup	2015
CERAFIT HAC cup	2004
Biomaterials	
CerafixGenta HV/BV cement	1985 (BV) et 2002 (HV)
Cerafix HV	2002

11. LOGOS

	Manufacturer		Date of manufacture		Medical device
	Batch code		Keep away from sunlight		CMR substance
	Consult instructions for use		Catalogue number of the product		Fears humidity
	IUD Code		Non-sterile		
	CE marking followed by the number of the notified body, (0459, LNE/G-MED)		Please note that the special precautions for use should be checked for any important safety information such as warnings or precautions to be taken		

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