Français

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Matériel ancillaire Rachis

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● INSTRUCTIONS FOR USE

Surgical instruments for the Spine Recommendations for reprocessing reusable

1. PURPOSE

These instructions are recommended for the upkeep, cleaning, maintenance and sterilization of manual and reusable orthopaedic surgical instruments for the CERADISC device manufactured and/or distributed by Ceraver®. This document has been created to assist healthcare staff in the safe handling, effective maintenance and reprocessing of Ceraver® reusable instruments.

2. SCOPE

This instruction manual provides information on the upkeep, cleaning, disinfection, maintenance and sterilisation of reusable surgical instruments for the CERADISC device manufactured and/or distributed by Ceraver®.

3. GLOSSARY

- **Contaminated:** state of an item having actually or potentially been in contact with micro-organisms.
- **Decontamination:** use of physical or chemical means to remove, deactivate or destroy blood-borne pathogens on a surface or item to the point where they can no longer transmit infectious particles and the surface or items are safe for handling or destruction.
- **Disinfection:** process used to reduce the number of viable micro-organisms on a product to a level previously specified as appropriate for handling or use.
- **Washer/Disinfector:** machine designed to clean and disinfect medical instruments and other items used in medical, dental, pharmaceutical, and veterinary practice settings.
- **Cleaning:** removal of any contamination on an item, required before any other processing.

Note: cleaning and disinfection are often carried out in a single operation.

- **Manual cleaning:** cleaning without using an automatic washer or washer/disinfector.
- **Sterile:** free from all viable micro-organisms.
- **Sterilisation:** validated process used to remove all forms of viable micro-organisms from an instrument.

Note: In a sterilisation process, the nature of microbiological death is described by an exponential function. Therefore, the presence of micro-organisms on an item can be expressed in terms of probability. Although this probability can be reduced to a very low number, it can never be reduced to zero. This probability can only be assured for validated processes.

- **Processing/Reprocessing:** activity consisting of cleaning, disinfection, and sterilisation, required to prepare a new or used surgical instrument prior to its use.

4. ACRONYMS

BI	= Biological Indicator
PPE	= Personal Protective Equipment
TSE	= Transmissible Spongiform Encephalopathy
CJD	= Creutzfeldt-Jakob Disease
SAL	= Sterility Assurance Level
SSD	= Sterile Supply Department
OT	= Operating Theatre
AORN	= Association of periOperative Registered Nurses

5. CONSIDERATIONS

This instruction manual refers to all Ceraver® reusable surgical instruments for the CERADISC device and must be studied carefully.

Ceraver® surgical instruments have been designed for implant placement and have been manufactured from materials whose biocompatibility is recognised for the use indicated.

The user should refer to and follow the surgical technique of the implant in question.

The user/controller must comply with local laws and regulations in countries where the reprocessing requirements are more stringent than those detailed in this manual.

New and used instruments must be thoroughly processed according to these instructions prior to use.

During musculoskeletal surgery, instruments are contaminated with blood, bone and tissue debris. They can also be contaminated with body fluids containing the hepatitis virus, HIV or other aetiological agents and pathogens. All healthcare personnel should be familiar with the universal precautions necessary to prevent injuries caused by handling instruments with cutting edges and/or sharp points during and after surgical procedures or during reprocessing.

It should be noted that physiological saline and other irrigation fluids are often used in large quantities during surgical procedures and exert a corrosive effect on the instruments.

Orthopaedic surgery requires heavy instruments comprising several components, articulated or rotating parts and removable handles. Instruments are often supplied in sets and divided into trays in which they are arranged by size or in the order required for a specific surgical procedure.

It is the responsibility of healthcare facilities to clean, disinfect, package and sterilise all instrument sets before returning them to Ceraver®. However, the next user must also inspect the set upon receipt to verify that the instruments have been properly cleaned and decontaminated before repeating the reprocessing procedures to prepare the loan set for subsequent reuse.

Core orthopaedic instrument sets must be complete and in good condition to be used correctly. Optional instruments may be available on request from your Ceraver® representative. To maintain the instruments correctly, it is important to consider the following information and processing instructions:

- Warnings and precautions
- Instrument set completeness and functionality
- Reprocessing limitations and/or restrictions
- Preparation for processing at the point of use

- Preparation for cleaning (including assembly/disassembly, if necessary)
- Cleaning, disinfection and drying
- Maintenance, inspection, testing and lubrication
- Sterile packaging
- Sterilisation
- Storage

6. PROCESSING INSTRUCTIONS

The surgical instruments are supplied non-sterile and must be cleaned and sterilised prior to use. After use, these instruments must, as a minimum, be decontaminated, cleaned, put away, and sterilised. The following information describes the correct steps, based on Ceraver®'s experience and the recommendations of organisations (see references), for the processing of Ceraver® surgical instruments to ensure they have a safe and efficient service life.

Note: these instructions describe the necessary processing steps that new and used instruments must undergo to attain sterility.

Warnings and precautions:

- **Universal precautions must be observed** by hospital personnel who work with contaminated or potentially contaminated medical instruments. Instruments with sharp points or cutting edges should be handled with extreme caution.
- **Personal Protective Equipment (PPE) must be worn** when handling or working with contaminated or potentially contaminated equipment, instruments and materials. PPE includes lab coat, mask, protective goggles or face shield, gloves and shoe covers.
- **Do not use metal brushes or scouring pads** during manual cleaning procedures. These materials may damage the surface and finish of the instruments. Use soft-bristled, nylon brushes and pipe cleaners.
- Cleaning agents with low foaming surfactants should be used during manual cleaning procedures to ensure that instruments are visible in the cleaning solution. Manual scrubbing with brushes should always be performed with the instrument immersed in the cleaning solution to prevent generation of aerosols and splashing, which may spread contaminants. Cleaning agents must be gently and thoroughly rinsed from instrument surfaces to prevent accumulation of detergent residues.
- **Do not place heavy instruments on top of delicate instruments.**
- **Do not allow contaminated instruments to dry before reprocessing.** All subsequent cleaning and sterilisation steps are facilitated by not allowing blood, body fluids, bone and tissue debris, saline solution or disinfectant to dry on used instruments.
- Saline solution and cleaning/disinfection agents containing aldehyde, mercury, active chlorine, chloride, bromine, bromide, iodine or iodide are corrosive and **must not** be used. Instruments **must not** be placed or immersed in **Ringer's solution**.
- Only instruments manufactured and/or distributed by Ceraver® should be included in Ceraver® instrument trays. These validated reprocessing instructions are not applicable to Ceraver® trays that include instruments that are not manufactured and/or distributed by Ceraver®.
- Descaling agents that include morpholine should not be used in steam sterilisers. These agents leave residues that can damage polymer instruments.
- It is important to take special care when manually cleaning instruments for boring, reaming or drilling so that all bone residue is properly removed.
- If an instrument breaks in situ, a fragment can be located using imaging techniques tailored to the material.

Inspection on receipt:

- Upon receipt, instrument sets should be inspected to check for completeness. In addition, illustrated bills of materials are provided with the trays to check the presence of instruments on arrival and departure from the facility.
- Surgical procedures follow a specific order in which the instruments are used. In addition, many instruments have dimensional features that aid in determining prosthesis size and height, etc. therefore, it is very important that all sizes required for a specific set of instruments are available. Contact your Ceraver® representative if requested instruments have been omitted but are required for the operation.

Limitations and restrictions:

Neutral pH enzymatic cleaning agents are recommended for cleaning Ceraver® reusable instruments. Alkaline agents with a pH of 12 or less can be used to clean stainless steel and some polymer instruments in countries where required by law or local regulations or where there is a risk of prion disease, such as Transmissible Spongiform Encephalopathy (TSE) and Creutzfeldt-Jakob Disease (CJD). **It is extremely important to neutralise and rinse alkaline cleaning agents from the instruments thoroughly and completely.**

Note: pins and other cutting instruments should be carefully inspected after processing with alkaline detergents to ensure that cutting edges are fit for use.

Note: It is important to select enzymatic solutions intended for breakdown of blood, tissues and body fluids. Some enzymatic solutions are designed specifically for breakdown of faecal matter or other organic contaminants and may not be suitable for use with instruments.

- Repeated processing in accordance with the instructions in this manual has minimal effect on Ceraver® reusable manual instruments, unless otherwise specified. End of life for stainless steel or metal surgical instruments is determined by wear and deterioration due to the intended surgical use and not to reprocessing.
- Automatic cleaning using a washer/disinfector alone may not be effective for orthopaedic instruments with lumens, cannulas, blind holes, contacting surfaces and other complex features. A manual/automatic cleaning process is recommended.
- Where applicable, multi-component instruments should be disassembled for cleaning. Disassembly, where necessary, is generally self-evident. Take care not to lose small parts. If a part is lost, notify your Ceraver® representative as soon as possible.
- Instruments must be removed from the metal or polymer trays for manual and/or automatic cleaning procedures. Do not clean instruments while they are in the polymer or metal trays. Instrument lids and trays must be cleaned separately from the instruments.
- The polymers used in Ceraver® instrument sets can be sterilised using steam/moist heat. Polymer materials have a limited shelf life (500 sterilisation cycles). If polymer surfaces become «chalky,» show excessive surface wear (e.g., crackling or delamination) or if the polymer instruments show excessive distortion or are visibly warped, they must be replaced. Notify the Ceraver® representative if polymer instruments need to be replaced.
- Most available polymers will not withstand conditions in washers/sterilisers that operate at temperatures of 141°C and above and use live steam jets for cleaning. Under these conditions, the polymer instruments will suffer severe surface damage.
- Immersion in disinfectants may be a necessary step to control certain viruses. However, these agents may discolour or corrode the instruments (domestic bleaches contain or form chlorine and chloride in the solution and have a corrosive effect similar to that of saline). Disinfectants containing glutaraldehyde or other aldehydes can denature protein based contamination agents, making them harder and more difficult to remove. Where possible, immersion in disinfectants should be avoided.
- Steam/moist heat is the recommended sterilisation method for Ceraver® instruments
- Ethylene oxide (EO) sterilisation, plasma sterilisation and dry heat sterilisation methods are not recommended for sterilisation of Ceraver® reusable instruments.
- Instruments with removable polymer sleeves **must** be disassembled for sterilisation.
- Avoid using acidic anti-corrosion agents.
- Use of hard water should be avoided. Softened tap water may be used for initial rinsing. Purified water should be used for final rinsing to eliminate mineral deposits on the instruments. One or more of the following processes may be used to purify water: ultrafiltration (UF), reverse osmosis (RO), deionisation or equivalent.

Point of use – preparation for processing:

- Remove excess body fluids and tissue from the instruments with a lint-free, disposable cloth. Place instruments in a basin of distilled water or in a tray covered with a damp cloth. Do not allow saline, blood, body fluids, bone and tissue fragments or other organic debris to dry on instruments prior to cleaning

Note: Soaking in proteolytic enzyme solutions facilitates cleaning, especially for complex instruments with hard-to-reach areas (e.g., tubular designs, etc.). These enzymatic solutions dissolve proteinaceous matter and prevent blood and protein-based materials from drying on the instruments. The manufacturer's instructions for preparation and use of these solutions must be followed explicitly.

- The instruments **must be** cleaned within 30 minutes of use to minimise the possibility of drying prior to cleaning.
- Used instruments **must be** transported to the supply department in closed or covered containers to prevent any unnecessary contamination risk.

Preparation before cleaning:

- Specific symbols or instructions engraved on instruments or trays must be strictly followed.
- Where applicable, multi-component instruments should be disassembled for appropriate cleaning. Care should be taken to avoid losing small screws and components. If a part is lost, notify your Ceraver® representative as soon as possible.
- Published instructions for use and surgical techniques and/or procedures may serve as an additional source of information to illustrate the assembly/disassembly instructions for specific Ceraver® instruments.

Preparation of cleaning agents:

- Ceraver® recommends using neutral pH enzymatic cleaning agents with low foaming surfactants. Alkaline agents with a pH of 12 or less can be used in countries where the laws or local regulations require it. Use of alkaline agents should be followed by a neutralising agent and thorough rinsing.
- All cleaning agents should be prepared in accordance with the dilution and temperature recommendations of the manufacturer. Softened tap water may be used to prepare cleaning agents. Use of the recommended temperatures is important for optimal performance of cleaning agents.
- Dry powder cleaning agents should be completely dissolved prior to use to prevent oxidation or corrosion of the instruments.
- Fresh cleaning solutions should be prepared when existing solutions become visibly contaminated (bloody and/or cloudy).

CLEANING PRODUCTS	
Water	Use cold demineralised water or reverse osmosis filtered water. Temperatures higher than 60°C (140°F) should be avoided, as they have a coagulating effect on proteins, making them difficult to remove.
Detergent	Prepare a detergent (e.g., CIDEZYME, LIQUI-NOX®, Alconox, Inc. pH 8.5) according to the manufacturer's instructions.
Enzymatic cleaner	Prepare an enzymatic cleaner (e.g., CIDEZYME, ENDOZYME®, Ruhof Corporation pH 6.0 to 7.5) according to the manufacturer's instructions.
Accessories for manual cleaning	Bristle brushes and/or bottle brushes, syringes, gloves, disposable absorbent cloths (e.g., KIMWIPE®, Kimtech Science)
Ultrasonic cleaner	Ultrasonic cleaners must be checked regularly to ensure that they are functioning correctly

Instructions for disinfection and manual pre-cleaning instructions before an automatic cycle

CLEANING/DISINFECTION		
Pre-Cleaning	Step 1	Immerse the instruments completely in an enzymatic solution and allow to soak for 10 minutes. Use a soft nylon bristle brush to gently brush the instrument until all visible dirt is removed. Pay particular attention to roughness, lumens, contact surfaces, connectors and other areas that are difficult to clean. Lumens should be cleaned with a long, narrow, soft nylon bristle brush.
	Step 2	Remove the instruments from the enzymatic solution and rinse them in purified water for at least 1 minute. Rinse lumens, holes and all other difficult-to-reach areas carefully and thoroughly.
	Step 3	Place instruments in the rack of a suitable washer/disinfector and process through a standard instrument washer/disinfector cycle.

Instructions for automatic disinfection/cleaning

Cleaning/disinfection	
Automatic cleaning	1. An automatic cleaning procedure can involve a washer-steriliser, a washer-disinfector, an ultrasonic washer, or other types of machines for washing and decontaminating instruments. There are many types of automatic washer systems, each with their own instructions that must be followed. 2. Washer-disinfector systems are not recommended as the sole cleaning method for surgical instruments. Instruments must be cleaned according to the combination of manual/automatic cleaning indicated in this manual, unless otherwise specified. 3. Reusable surgical instruments for the Ceraver® CERADISC can be cleaned and disinfected using a typical washer/disinfector cycle including at least: - pre-cleaning <45° C - 2 minutes - cleaning at 55°C - 5 minutes - neutralization/rinse phase - 2'2" minutes - a thermal disinfection phase at 90°C - 5 minutes At the end of this cycle, instruments should be carefully inspected.

Inspection, maintenance, testing and lubrication

1. Carefully inspect each instrument to ensure that all visible contamination has been removed. If contamination is found, repeat the cleaning/disinfection process.
2. Visually inspect for completeness, absence of damage and/or excessive wear.

Note: If there is damage or wear that could compromise the function of the instrument, contact the Ceraver® representative for a replacement.

3. Verify proper functioning of moving parts (e.g., safety tip, nut, sliding parts, etc.).
4. Hinged, rotating or articulated instruments should be lubricated with a water soluble product (e.g., instrument milk or an equivalent lubricant) intended for surgical instruments to be sterilised.

Note: Do not use mineral oil or silicone lubricants, as they:

- a. coat micro-organisms;
- b. prevent direct contact of steam with surfaces;
- c. are difficult to remove.

Note: These instructions do not apply to electric or pneumatic instruments. These instruments have different requirements and should be lubricated according to the manufacturer's instructions.

5. Check for distortion of elongated shaped instruments.
6. Optional Ceraver® instruments should not be added to a preconfigured instrument tray unless a dedicated universal space or compartment has been included in the design and the guidelines described below for trays without defined layouts or universal spaces can be applied.

Packaging:

Packaging of instrument sets in rigid trays with lids.

Instrument trays with defined, preconfigured layouts.

- Areas designated for specific instruments must contain only instruments specifically intended for these areas.
- Optional Ceraver® instruments should not be added to a preconfigured instrument tray unless a dedicated universal space or compartment has been included in the design and the guidelines described below for trays without defined layouts or universal spaces can be applied.
- Only instruments manufactured and/or distributed by Ceraver® should be included in Ceraver® instrument trays. These validated reprocessing instructions are **not applicable** to Ceraver® trays that include instruments that are not manufactured and/or distributed by Ceraver®.

Universal instrument cases and trays without defined, preconfigured layouts or containing undefined universal spaces or compartments should only be used under the following conditions:

- If possible, all instruments must be disassembled before being placed in the tray.
- Instructions for assembly and disassembly are included with the corresponding bill of materials and sent to the operating theatre for assembly/disassembly of disassembled sterilised instruments.
- All the instruments should be arranged to ensure steam penetration to all instrument surfaces. Instruments should not be stacked or placed in close contact.
- The user must check that the instrument tray has not tipped or the contents shifted once the instruments have been arranged in the tray.
- Only instruments manufactured and/or distributed by Ceraver® should be included in Ceraver® instrument trays. Ceraver® validated reprocessing instructions **are not applicable** to Ceraver® trays that include instruments that are not manufactured and/or distributed by Ceraver®.

Sterilisation instructions:

- Refer to the table below for the recommended minimum sterilisation parameters that have been validated by Ceraver® to offer a Sterility Assurance Level (SAL) of 10⁻⁶.

Cycle type	Minimum temperature	Minimum exposure time	Minimum dry time
Prevacuum	134°C	18 min	30 min

Note: the steriliser manufacturer's instructions for operation and load configuration should be followed explicitly.

- The facility is responsible for the internal procedures for the reassembly, inspection and packaging of instruments once they have been thoroughly cleaned in a manner that will ensure sterilising steam penetration and adequate drying. Provisions for protection of any sharp or potentially dangerous edges of the instruments should also be recommended by the care facility.
- Steam/moist heat sterilisation is the recommended method for Ceraver® instrument sets.
- Instrument sets should be properly prepared and packaged in Ceraver® trays designed to allow steam to penetrate and make direct contact with all surfaces.
- Ethylene oxide or plasma sterilisation methods **should not** be used.
- Gravity displacement sterilisation cycles are not recommended because cycle times are too long.

Storage instructions:

- Sterile, packaged instruments should be stored in a designated limited access area that is well ventilated and protected from dust, moisture, insects, and extreme temperatures/humidity.
- Sterile instrument packages should be carefully examined prior to opening to ensure that package integrity has not been compromised.

Note: maintenance of sterile package integrity is generally event related. If sterile packaging is torn, perforated, shows evidence of tampering or has been exposed to moisture, the instrument set must be repackaged and sterilised.

Note: if there is any evidence that the lid seal or filters on a sterilisation container have been opened or handled, the sterile filters must be replaced and the instrument set re-sterilised.

7. FACILITIES' RESPONSIBILITIES FOR CERAVER® INSTRUMENTS LOAN SETS

- Orthopaedic surgical instruments generally have a long service life; however, mishandling or inadequate protection can rapidly reduce it. Instruments which no longer perform properly because of long use, mishandling or improper maintenance should be returned to Ceraver®. Notify your Ceraver® representative of any problem with the instruments.
- Loan sets should undergo all the steps of decontamination, cleaning, disinfection, inspection and final sterilisation before being returned to Ceraver®. Documentation of decontamination should be provided with instruments being returned to Ceraver®.
- In the event that instruments are found to be missing from or damaged in loan sets, the information should be reported to the OT contact person who will forward it to Ceraver® Sales Administration and your Ceraver® representative. This is necessary to ensure the traceability of the event, assess any costs generated and inform the Ceraver reception department to correct the instrument set and make it compliant for a future use.
- The instructions provided in this manual have been validated in the laboratory by Ceraver® and are fully tailored to the preparation of instruments for use. It is the responsibility of the facility to ensure that reprocessing is performed using the appropriate equipment and materials, and that the personnel tasked with the reprocessing have been correctly trained in order to achieve the desired result. The equipment and processes should be validated and monitored regularly. Any deviation from these instructions by the processing supervisor should be properly assessed for effectiveness to avoid any potential adverse consequences.

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9. REFERENCES

1. ISO 17664, Sterilisation of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices
2. ISO 17665-1, Sterilisation of Health Care Products – Moist heat

Compliance with standards ISO 17664 and ISO 17665-1 is noted in the sterility and cleaning validation reports.

10. NOTE

The CERADISC device ancillaries were CE marked in 2016.

11. LOGOS

	Manufacturer		Date of manufacture
	Batch code		Keep away from sunlight
	Consult instructions for use		Catalogue number of the product
	IUD Code		Non-sterile
	Medical device		Fears humidity
	CE marking followed by the number of the notified body, (0459, LNE/G-MED)		Please note that the special precautions for use should be checked for any important safety information such as warnings or precautions to be taken

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