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Français

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INSTRUCTIONS FOR USE

English

Instructions for use for the following implants:

- Bipolar heads (references: 2808, 2809, 2810, 2811, 2012, 2013, 2014, 2015, 2016)

Graphic symbols used for labelling these devices:

Symbol	Symbol title	Symbol	Symbol title
	Manufacturer		Do not reuse
	Expiry date		Product terminally sterilised by irradiation.
	Batch code		Indicates a medical product having undergone sterilisation.
	Catalogue number of the product		Caution, all the information required for the use of the device is not found on the label; it is necessary to consult the instructions for use.
	Do not use if the packaging is damaged		Do not resterilise
	CE marking followed by the code for the notified body (0459, LNE/G-MED)		Consult instructions for use
	Moisture sensitive		Keep away from sunlight

Before using a CERAVER® hip prosthesis, the attention of the surgeon is drawn to the following recommendations:

1- GENERAL POINTS

- The manufacturer and supplier of a CERAVER® hip prosthesis guarantee the quality of the manufacture and materials. A CERAVER® hip prosthesis that has already been used must not be re-implanted. Instructions for the assembly of the different components must be adhered to as defined in the commercial and technical operating documentation. The use of CERAVER® prosthetic hip components with components of other origins is strictly prohibited.
- Several criteria are of crucial importance for the success of an operation using a CERAVER hip prosthesis:
 - Appropriate patient selection for the relevant procedure.
 - Perfect operating technique and correct use of the surgical instruments specially designed by CERAVER®.
 - Application of strict aseptic technique.
- The surgeon is responsible for any complications that may arise from an erroneous indication and/or a defective operating technique or incorrect asepsis; these cannot be attributed to the manufacturer or the supplier of the CERAVER® products.
- The surgeon must be aware of all problems posed by implanting a joint prosthesis. The operating technique for implantation of a CERAVER® hip prosthesis may be taught in a surgical centre where surgeons have acquired perfect mastery of this technique.
- Patients receiving a prosthetic hip joint should be informed that their weight and activity level may affect the longevity of the prosthesis.

2- IMPLANT DESCRIPTION

CERAVER bipolar heads are manufactured using the latest technologies. The mechanical properties and functional properties can only be guaranteed if the surgical technique and present instructions are strictly observed.

- The bipolar heads consist of one part made of high nitrogen stainless steel (Standard ISO 5832-9) and two parts made of ultra-high molecular weight polyethylene (UHMWPE) (Standard ISO 5834-2).

3- INDICATIONS

- A hip prosthesis should be considered only after all other surgical possibilities have been carefully considered and deemed to be less appropriate.
- Even the most successfully implanted artificial joint is inferior to a sound natural joint. On the other hand, an artificial joint can be a highly beneficial substitute for a deformed or diseased joint, because it eliminates pain and is conducive to the restoration of good mobility and weight-bearing capacity of the operated limb.
- Every artificial joint is always subject to wear which cannot be avoided. An initially stable prosthesis may become loose in the course of time. Wear and loosening of a prosthesis may render a re-operation necessary.
- Infection of a prosthesis will usually necessitate its removal and any such infection may leave debilitating sequelae.

Approved applications:

Hip hemiarthroplasty involves partial replacement of a hip joint affected by various disorders, which cannot be treated by other therapies or methods. Its purpose is to give the patient back all or part of his/her mobility and stability while minimising or even eliminating the pain. Partial replacement of the hip joint may be indicated in cases of traumatic fracture of the head or neck of femur.

The surgeon must inform patients of the risks of implanting a prosthesis and their agreement to the suggested operation must be obtained.

The following circumstances require special precautions:

- Overweight patient.
- Excessive weight-bearing expected owing to profession or sports practices.
- Epilepsy or other factors that may cause repeated accidents presenting an increased risk of fractures.
- Severe osteoporosis and/or osteomalacia.
- Old infectious diseases with potential for evolution, with possibility of joint manifestations.
- Reduction in bone support capacity as a result of a tumour.
- Intellectual incapacity of patients to understand and follow the surgeon's instructions.
- Alcoholism or drug addiction;
- Metabolic disorders or general drug treatments leading to progressive deterioration of the bone support for the prosthesis;
- Disorders of other joints (knee, ankle, etc.);
- Overall weak immune condition of the patient.

4- CONTRAINDICATIONS

- The following circumstances generally constitute an absolute contraindication to the implantation of a joint prosthesis:
- Allergy to any component of the prostheses.
 - Patients with a native acetabulum that has a diameter of less than 40 mm or greater than 56 mm.
 - Evolving infection of the joint or relevant periparticular area, including severe neuroarthropathy;
 - Loss of bone substances above or below the joint, making it impossible to anchor the cementless prosthesis correctly;
 - Muscle loss, neuromuscular lesion or vascular insufficiency of the affected limb, making the procedure unjustifiable;
 - Charcot's disease;
 - Patients under 50 years old with joint lesions that could probably be treated successfully by a conservative treatment;
 - Pregnancy;
 - Tumours.

5- IMPORTANT WARNINGS

- 5.1. Preoperative warnings**
 - The operation must be scheduled with care. Depending on the results of the X-ray examination (templates) and preoperative planning, the surgeon will decide on the suitable implant.
 - CERAVER® instruments and templates specially designed to prepare for the implantation of the joint components should be used. The description of the surgical procedure can never be exhaustive or specify all the risks and complications that must be accounted for. Technical notes on our products and publications with more detailed information, as well as video documents, may be requested from CERAVER® or its distributor.
- 5.2. Perioperative warnings**
 - The joint prosthesis must be handled with care. Scratches or unevenness due to impacts will lead to excessive wear and can result in an increased risk of complications.
 - For the choice of prosthesis size: CERAVER templates, ancillary instruments and test implants will be the only ones used.
 - Exposed prosthetic joint surfaces must not be in contact with any metal object, instruments for example.
 - The insert and supporting head must be assembled in accordance with the surgical procedure and using CERAVER instruments.
 - The various components of the CERAVER bipolar head must be used together and require prior assembly with a CERAVER femoral head before connecting to the CERAVER femoral stem.

6- POSTOPERATIVE PRECAUTIONS

- All patients implanted with a CERAVER® joint prosthesis should see an orthopaedic surgeon for postoperative follow-up. This examination is important to detect any signs of prosthetic wear or loosening before the appearance of functional manifestations. Early recognition of a complication enables effective measures to be taken. In particular, early re-operation has more chance of success than a late operation, particularly when advanced loosening has already led to significant resorption lesions in the bone.
- The patient must inform his/her surgeon about the slightest change occurring in the operated limb.
- If an annual follow-up appointment is not possible, a control X-ray image should be sent to the surgeon for evaluation at least every year. The radiological study makes it possible to screen for possible early lesions.
- 6.4. MRI compatibility.** The components of CERAVER® bipolar heads are non-ferromagnetic. They have a geometry that is unlikely to generate eddy currents, and are fixed to the bone or the implant, reducing the risk of mobilisation. Their safety, particularly in terms of implant heating and migration, was assessed from bibliographic data compared with available data on devices of similar composition, shape and uses. This evaluation has concluded that there is no danger for MRI examinations. This evaluation has concluded there is conditional compatibility for MRI examinations up to 3 Tesla. However, as a precaution, it is recommended that MRI examinations are avoided in the 48 hours after implantation, and that the person responsible for the examination of the recent insertion of an implant is warned, if such an examination proves to be essential. It should be noted that all devices presenting a high contrast with the surrounding biological environment may generate artefacts which must be taken into account for the perfect formation and interpretation of imaging examinations. For this purpose, a patient with an implant should be advised to inform the health professionals concerned (radiologists and radiology technicians) about the presence of an implant as soon as possible before these examinations.

7- COMPLICATIONS

The following complications have occasionally been observed, and require special attention:

- Loosening of the prosthesis following changes in weight-bearing transmission or tissue reactions due to the prosthesis. Loosening is frequently the result of one or more of the risk factors previously mentioned.
- Early and late-onset infections.
- Dislocation, subluxation, insufficient range of motion, undesirable shortening or extension of the operated limb following a non-optimal implantation of the prosthesis.
- Bone fractures following uniaxial overload or bone fragility.
- Transient or permanent nerve lesions following compression or haematoma.
- Cardiovascular diseases, including vein thrombosis, pulmonary embolisms and cardiac arrests.
- Haematoma of the surgical wound and delay in healing.
- Tissue reactions following an allergy to the implanted materials, particularly metal, or following the accumulation of wear particles.
- Temporary or permanent peripheral neuropathy;
- Pulmonary disorders, including pneumonia and atelectasis;
- Urological complications (retention, urinary infection);
- Other associated complications resulting from anaesthesia, use of ancillary equipment, nosocomial diseases;
- Ectopic ossification;
- Pain.

8- COMPATIBLE MEDICAL DEVICES

- The stainless steel bipolar heads are compatible with:
 - Dedicated ancillary equipment manufactured by CERAVER;
 - The following implants manufactured by CERAVER:
- Stainless steel femoral heads and alumina femoral heads are compatible with bipolar heads of the same joint diameter.
- It should be noted that bipolar heads with reference 2808 to 2811 are compatible only with stainless steel heads of 22.2 mm in diameter, and bipolar heads with reference 2012 to 2016 are compatible only with heads of 28 mm in diameter.

9- STORAGE AND HANDLING PROSTHESES

- Joint prostheses must be stored in the original packaging.
- Joint prostheses must not be subjected to treatment or modification.
- Prostheses that are damaged or have undergone a treatment that is inappropriate or not authorised by the manufacturer must not be implanted.
- Do not use the device after the expiry date. Contact your CERAVER® representative and return the product to him/her.

10- PACKAGING AND STERILISATION

The medical devices are packed individually in double packaging and are supplied sterile. CERAVER® bipolar heads are sterilised by gamma radiation at a minimum dose of 25 kGy. It is important to check the integrity of the sterility protect before use.

11- RESTERILISATION AND REUSE

Resterilisation and/or reuse of any of the components of a CERAVER® hip prosthesis is strictly prohibited at the risk of exposing the patient concerned to contamination that could cause an infection. Reuse of the constituent parts of a CERAVER hip prosthesis would lead to a degradation in performance of the devices concerned (poorer mechanical resistance, increased risk of implant fracture).

12- DISPOSAL AFTER USE

If a device is removed for revision surgery, dispose of the components according to the centre's medical waste disposal procedure. This does not apply if the follow-up procedure is due to a Medical Device Vigilance case.

13- NOTE

Bipolar heads were CE marked in 1995. Federal law in the USA and Canada requires that prosthetic components be purchased by a pharmacist. For further information concerning the use of these components, please contact the manufacturer, their overseas subsidiaries and distributors.

CERAVER - Les Laboratoires Osteal Médical
69 rue de la Belle Etoile
95957 Roissy CDG Cedex - France
Tel. : 33 (0)1 48 63 88 63 - Fax : 33 (0)1 48 63 88 99
www.ceraver.com

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