Français

English

● NOTICE D'INSTRUCTIONS

SUBSTITUTS OSSEUX CALCIRESORB

● INSTRUCTIONS FOR USE

CALCIRESORB BONE SUBSTITUTES

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CALCIRESORB-KNOCHENERSATZ

● MANUAL DE INSTRUCCIONES

SUSTITUTOS ÓSEOS CALCIRESORB

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SOSTITUTI OSSEI CALCIRESORB

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SUBSTITUTS OSSEUX CALCIRESORB

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INSTRUCTIONS FOR USE

CALCIRESORB BONE SUBSTITUTES

CALCIRESORB *Classic*[®]
CALCIRESORB *Bone Like*[®]
CALCIRESORB *Collagen*[®]

Document pertaining to the following device(s):

- CALCIRESORB *Classic*[®]** bone substitutes (Ref. 1612-1613-1614-1615-1616-1627-1628-1617-6620-6621-1604-1605-1606-1607-1608-1610-6700-6701-6702-6705-6706-6707-6708-6730-6731-6734-6735-6736-6740-6741-6741-6742-6743-6744)
- CALCIRESORB *Bone Like*[®]** bone substitutes (Ref. 6930-6931-6932-6933-6934-6935-6937-6938-6945-6946-6947-6948-6949)
- CALCIRESORB *Collagen*[®]** bone substitutes (Ref. 1623-1629-1630-1631)

Recommendations

Before using a **CALCIRESORB[®]** bone substitute, surgeons must familiarise themselves with the following recommendations:

1/ GENERAL POINTS

1.1. The manufacturer and supplier of a **CALCIRESORB[®]** bone substitute guarantee the quality of the manufacturing processes and materials used.

1.2. Several criteria are key to the success of a procedure using a calcium-phosphate bone substitute:

- careful selection of a suitable patient for the operation,
- extremely strict asepsis.

1.3. The surgeon is responsible for any complications that may arise due to an erroneous indication, a faulty operating technique or poor asepsis. Neither the manufacturer nor the supplier of the CERAVER[®] products may be held liable for any such complications.

2/ MANUFACTURING – MATERIALS

Implants in the **CALCIRESORB *Classic*[®]** and **CALCIRESORB *Bone Like*[®]** ranges are manufactured out of tricalcium phosphate β (β - $\text{Ca}_3(\text{PO}_4)_2$; synthetic origin), using the most recent technologies to ensure excellent purity quality. Their mechanical and functional properties can only be guaranteed if the present instructions are strictly observed.

The **CALCIRESORB *Collagen*[®]** composites are produced using **CALCIRESORB *Classic*[®]** beads coated with a matrix of highly purified bovine collagen fibres. The combination of the chemical composition (pure TCP) and the specific structural properties enables the colonisation of macropores by the newly formed bone, together with gradual but total resorption of the **CALCIRESORB[®]** implants:

CALCIRESORB *Classic*[®] / CALCIRESORB *Collagen*[®] :

- 30% porous materials ($\pm$ 5%): the size of the macropores ranges from 100 μm to 400 μm , giving an interconnection rate of over 80%.
- 45% porous materials ($\pm$ 5%): the size of the macropores ranges from 100 μm to 400 μm , giving an interconnection rate of over 80%.

- CALCIRESORB *Bone Like*[®] :**
- 75% porous materials ($\pm$ 5%): the macropores have an average size of 500 μm , giving an interconnection rate of over 96%.

3/ STORAGE AND HANDLING

- The bone substitutes must be stored in their original packaging.
- The bone substitutes must not be modified in any way, except where authorised by these instructions.
- Damaged bone substitutes must not be implanted.
- Do not exceed the use-by date indicated on the packaging.

4/ PACKAGING AND STERILISATION

- The bone substitutes are delivered in sterile double packaging.
- They are sterilised by gamma radiation at a minimum dose of 25 kGy;
- It is important to check the integrity of the sterility protection before use.

5/ RESTERILISATION AND REUSE

The resterilisation and/or reuse of the **CALCIRESORB[®]** bone substitutes (***Classic*[®]**, ***Bone Like*[®]** or ***Collagen*[®]**) is strictly prohibited, as it could expose the patient concerned to contamination and a subsequent risk of infection. The resterilisation and/or reuse of a **CALCIRESORB[®]** bone substitute would reduce its efficiency, notably through reduced mechanical strength (risk of fracture) and an incomplete bone remodelling phenomenon. There is also a risk of severe inflammatory reaction (for example, a reaction against a pre-implanted graft).

6/ INDICATIONS

CERAVER bone substitutes are used as spongy or cortico-spongy bones:

- for post-traumatic reconstruction,
- for bone-void filling (loss or insufficiency of substance),
- to assist fusion in arthrodesis (vertebral, talocrural, radiocarpal, etc.).

The surgeon chooses which implant to use. Provided that CERAVER's recommendation to avoid putting the products under mechanical stress without osteosynthesis equipment is observed, the applications described below are feasible using implants from the following ranges: **CALCIRESORB *Classic*[®]**, **CALCIRESORB *Bone Like*[®]**, and **CALCIRESORB *Collagen*[®]**, whether in the form of parts or beads:

- Orthopaedic and Traumatological Surgery**: bone-void filling during the implantation or revision of osteoarticular prostheses; assisting with the consolidation of fractures and pseudarthroses; articular and vertebral arthrodeses (intersomatic, posterolateral or intra-articular) in degenerative spine disease or trauma; bone-void filling after tumoral resection and curetting; osteotomies; restoration of the tibial plateau; and restoration of the calcaneum.
- Maxillofacial and Dental Surgery**: filling to assist consolidation following facial trauma; elevation of the sinusal floor; post-extraction alveolar filling; and implant reinforcement.
- Neurosurgery**: filling of trepanation holes; and spinal surgery (see "Orthopaedic and Traumatological Surgery" above).

The use of beads and composites is recommended for filling irregularly shaped cavities. The use of **CALCIRESORB *Collagen*[®]** composites is indicated in the case of cavities that cannot be closed, to prevent the beads from escaping the cavity.

The **CALCIRESORB[®]** range is presented below.

CALCIRESORB *Classic*[®]

- Implant for the restoration of tibial plateaux: corners.
- Implant for tibial or femoral osteotomies: oval-shaped corners.
- Implant for void filling: parallelepipeds; cylinders.
- Implant for arthrodesis: parallelepipeds.
- Spheroidal beads for void filling: 4 sizes available (0.5 - 1 mm, 1 - 2 mm, 2 - 3 mm, 3 - 4 mm) for volumes ranging from 2 to 30 cm³.

CALCIRESORB *Bone Like*[®]

- Implant for arthrodesis or void filling: parallelepipeds; cylinders.
- Cubic beads for void filling: 1 size available (4 mm) for volumes ranging from 2 to 30 cm³.

CALCIRESORB *Collagen*[®]

- Implant for arthrodesis or void filling: parallelepipeds; cylinders.

7/ RECOMMENDATIONS

- The surgeon must choose the graft size and shape best suited to the surgical procedure concerned. Radiological examination will facilitate the surgeon's choice.
- The osseous site that will be in contact with the graft must be cleaned and, if possible, well vascularised to ensure that a high-quality bone stock is available.
- For the graft to be fully effective, it must be placed in close, stable contact with the osseous site (by setting system or osteosynthesis). Additionally, and to improve the graft, it is highly advisable to impregnate the implant with bone marrow or, failing that, with autologous blood, and to add, where possible, spongy chips from drilling, grinding or cleaning.
- The bone substitutes should only be used in areas of mechanical stress when combined with osteosynthesis equipment capable of ensuring the stability of the surgical assembly, as is the case for any reconstruction with bone of autogenous or allogeneous origin.
- The bone substitutes in the **CALCIRESORB[®]** range are compatible with an MRI examination with no artefact risk.

8/ CONTRAINDICATIONS

The following circumstances constitute a contraindication to bone grafting:

- general contraindications to grafting;
- acute osteomyelitis or necrosis of the host site;
- bone degeneration;
- known active infection;
- allergy to any constituent of the bone substitute;
- cases where any surgery is contraindicated.
- Under no circumstances should **CALCIRESORB *Collagen*[®]** bone substitutes be implanted in direct contact with the central nervous system (brain, spinal cord).

9/ EXPECTED PERFORMANCE

As their chemical composition is very close to that of natural mineral bone, calcium phosphates – and particularly β -tricalcium phosphate – play an important role in the biological processes governing the stability and regeneration of bone tissue. Accordingly, the role of **CALCIRESORB[®]** bone substitutes is to temporarily replace bone tissue in the event of major tissue deficiency caused by trauma or a physiological problem. Thus these devices enable initiation of the bone remodelling phenomenon whilst facilitating their replacement with newly formed natural bone.

10/ DATE OF FIRST EC MARKING

- CALCIRESORB *Classic*[®]** bone substitutes: 1996
- CALCIRESORB *Bone Like*[®]** bone substitutes: 2007
- CALCIRESORB *Collagen*[®]** bone substitutes: 1996

11/ DETAILS FOR THE USE OF DEVICES IN THE UNITED STATES AND CANADA

In the USA and Canada, federal law requires that prosthetic components be purchased by a pharmacist. For further information concerning the use of these components, please contact the manufacturer or its foreign subsidiaries and distributors.

DATE OF LAST REVISION OF THIS DOCUMENT: **September 2013**

CERAVER - Les Laboratoires Osteal Medical
69, rue de la Belle Etoile – CS 54263
95957 ROISSY CDG cedex
Tél : +33 1 48 63 88 63
Fax : +33 1 48 63 88 99
www.ceraver.com

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