

**COMPOSANT FEMORAL HiFit PS**


Taille	AP (mm)		ML	Gauche	Droit
	Ext	Box			
1	51,5	31,6	60	15501	15601
2	55	34,1	65	15502	15602
3	57,5	36,6	65	15503	15603
4	60	39,1	70	15504	15604
5	62,5	41,6	70	15505	15605
6	65	44,1	75	15506	15606
7	67,5	46,6	75	15507	15607
8	70	49,1	80	15508	15608
9	75	54,1	85	15509	15609

**PLATEAU TIBIAL HiFit PS**


PLOT	TAILLE	AP (mm)	ML (mm)	REFERENCE					
				EPAISSEUR (en mm)					
				7	9	11	13	15*	17*
PETIT PLOT (ø 17mm)	A	39	63	15907	15909	15911	15913	15915	15917
	B	42	66,5	15927	15929	15931	15933	15935	15937
	C	45	70	15947	15949	15951	15953	15955	15957
	D	48	74	15987	15989	15991	15993	15995	15997
GRAND PLOT (ø 21,7mm)	C	45	70	15967	15969	15971	15973	15975	15977
	D	48	74	16007	16009	16011	16013	16015	16017
	E	52	78,5	16027	16029	16031	16033	16035	16037
	F	56	84	16047	16049	16051	16053	16055	16057

**EMBASE TIBIALE HiFit PS**


Taille	AP (mm)	ML (mm)	REFERENCE
<b>A</b>	39	63	<b>15701</b>
<b>B</b>	42	66,5	<b>15703</b>
<b>C</b>	45	70	<b>15705</b>
<b>D</b>	48	74	<b>15707</b>
<b>E</b>	52	78,5	<b>15709</b>
<b>F</b>	56	84	<b>15711</b>

**BOUTON ROTULIEN**


Ø en mm	épaisseur en mm	REFERENCE
<b>31</b>	8	<b>5090</b>
<b>34</b>	9	<b>5091</b>
<b>36</b>	10	<b>5092</b>
<b>39</b>	11	<b>5093</b>

**RALLONGE TIBIALE**


Ø DE LA RALLONGE en mm	LONGUEUR DE LA RALLONGE en mm	REFERENCE
<b>10</b>	50	<b>5080</b>
<b>11</b>	50	<b>5081</b>
<b>12</b>	50	<b>5082</b>
<b>13</b>	50	<b>5083</b>
<b>14</b>	50	<b>5084</b>
<b>15</b>	50	<b>5085</b>
<b>16</b>	50	<b>5086</b>
<b>17</b>	50	<b>5087</b>
<b>18</b>	50	<b>5088</b>
<b>19</b>	50	<b>5089</b>


**REFERENCE**
**3051**