

CEMENTED PS FEMORAL COMPONENT - HiFit PS


SIZE	AP (mm)		ML(mm)	LEFT	RIGHT
	Ext	Box			
1	51,5	31,6	60	15501	15601
2	55	34,1	65	15502	15602
3	57,5	36,6	65	15503	15603
4	60	39,1	70	15504	15604
5	62,5	41,6	70	15505	15605
6	65	44,1	75	15506	15606
7	67,5	46,6	75	15507	15607
8	70	49,1	80	15508	15608
9	75	54,1	85	15509	15609

PS TIBIAL INSERT - HiFit PS


PEG	SIZE	AP (mm)	ML (mm)	REFERENCE					
				TICKNESS (in mm)					
				7	9	11	13	15*	17*
SMALL PEG (\varnothing 17mm)	A	39	63	15907	15909	15911	15913	15915	15917
	B	42	66,5	15927	15929	15931	15933	15935	15937
	C	45	70	15947	15949	15951	15953	15955	15957
	D	48	74	15987	15989	15991	15993	15995	15997
BIG PEG (\varnothing 21,7mm)	C	45	70	15967	15969	15971	15973	15975	15977
	D	48	74	16007	16009	16011	16013	16015	16017
	E	52	78,5	16027	16029	16031	16033	16035	16037
	F	56	84	16047	16049	16051	16053	16055	16057

CEMENTED PS TIBIAL BASEPLATE - HiFit PS



SIZE	AP (mm)	ML (mm)	REFERENCE
A	39	63	15701
B	42	66,5	15703
C	45	70	15705
D	48	74	15707
E	52	78,5	15709
F	56	84	15711

CEMENTED MONOBLOC PATELLAR COMPONENT



Ø in mm	thickness in mm	REFERENCE
31	8	5090
34	9	5091
36	10	5092
39	11	5093

TIBIAL EXTENSION



Ø in mm	LENGTH in mm	REFERENCE
10	50	5080
11	50	5081
12	50	5082
13	50	5083
14	50	5084
15	50	5085
16	50	5086
17	50	5087
18	50	5088
19	50	5089

SCREW FOR TIBIAL EXTENSION



REFERENCE

3051